

Tracheostomy/ Laryngectomy

PRODUCT CATALOG



primed®



Primed[®] is a highly experienced manufacturer of medical devices for tracheostomized and laryngectomized patients.

For many years we have been supplying our customers with high quality products "MADE IN GERMANY".

In line with our corporate policy, the medical devices offered by **Primed**[®] are developed, produced and marketed in our own company. This is done in close cooperation with our customers and academic institutions.

The basis is formed by a certified development process in accordance with the DIN 13485 standard, in which highly innovative medical devices are continuously developed on the basis of ongoing market research and analysis.

Our vertical integration extends throughout the entire manufacturing process (injection moulding, extrusion, assembly) and also includes the sterilisation of devices in our own EO sterilisation unit.

The bundling of all the product-related processes in one company allows us to react flexibly to market requirements and guarantee the continuous further improvement of our existing products as well as to launch new products on the market.

PercuTrach®

PDT SETS

(Complete sets for percutaneous dilatational tracheostomy including tracheostomy tube) 7

 PDT Set Type 1 8

 PDT Set Type 2 14

Dilation sets (Sets for percutaneous dilatational tracheostomy without tracheostomy tube) 20

Accessories 21

Optima PercuTrach®

Tracheostomy tubes 22

Optima

Tracheostomy tubes 29

Replacement inner cannulas 38

Ultrasoft

Tracheostomy tubes 39

Replacement inner cannulas 48

Priflex®

Tracheostomy tubes 49

Priflex® MRI

Tracheostomy tubes 57

Primedistom®

Tracheostomy tubes cuffless 62

Tracheostomy tubes with Cuff 83

Replacement inner cannulas 90

Silver

Tracheostomy tubes 92

Accessories 109

NEO PED PEDL PEDC PEDCL

Tracheostomy tubes for Neonates and Children 110

 Neonatal tracheostomy tubes without Cuff 113

 Pediatric tracheostomy tubes without Cuff..... 114

 Pediatric tracheostomy tubes without Cuff, long..... 115

 Pediatric tracheostomy tubes with Cuff 116

 Pediatric tracheostomy tubes with Cuff, long 117

Accessories 118

Primedisilk®

Laryngectomy tubes 119

Tracheostomy placeholders 127

 Tracheosafe..... 128

 T-Tubus 129

 Primedibutton® 130

Primedistrip®

..... 133

Primed Tracheo Seal®

..... 139

Accessories for pulmonary rehabilitation 140

 Heat and Moisture Exchangers (HME) 141

 Accessories 147

 Speaking valves 148

 • for tracheostomized patients 149

 • for ventilated tracheostomized patients 153

 • for laryngectomized & tracheostomized patients 154

 • with filter function for laryngectomized & tracheostomized patients 155

Bacteria and virus filters 156

Tracheal compresses 157

Neck straps 158

Cough cap and shower accessories 160

Accessories for suctioning 161

Accessories for cleaning 162

Further accessories 163

Graphics 164

Pictograms 166

PDT SETS

(Complete sets for percutaneous dilatational tracheostomy including tracheostomy tube)

Percutaneous dilatational tracheostomy is a standardized method in intensive medical treatment. As alternative for surgical tracheostomy it may be carried out less time-consuming at the bedside in ICU's. The dilated tracheostoma closes spontaneously and long-term cosmetic results are far superior to conventional tracheostomy.

All **Primed®** PDT sets include a dilation set and a matching tracheostomy tube set .

PDT set type 1 includes a dilation set for the Ciaglia dilation technique and an Optima PercuTrach® tracheostomy tube with atraumatic inserter.

PDT set type 2 includes a dilation set for the Griggs dilation technique, with tracheal dilating forceps, and an Optima tracheostomy tube.

Our dilation sets are ideally adapted to the use of **Optima** and **Optima PercuTrach®** tracheostomy tubes which are available with or without sieving and with or without subglottic suction line.

Dilation sets and tracheostomy tube sets are also available separately.



DILATATION TRACHEOSTOMY

developed by Ciaglia


The percutaneous dilatation tracheostomy technique according to Ciaglia is a minimal-invasive treatment which implies the stepless bouginage of the trachea by a conically shaped dilator via guide wire, thus enabling the insertion of a suitable tracheostomy tube.

Primed® PDT Set Type 1 includes a conically shaped dilator which has a hydrophilic coating for assuring good gliding capabilities. After removing the dilator, the tracheostomy tube can be gently placed in the trachea via the atraumatic inserter.


The procedure should be performed under bronchoscopic control.




Position patient with hyperextension of neck, cover area around the neck with sterile dressing, disinfect skin




Palpation of neck, sonography (recommended), diaphanoscopy of neck area (obligatory), horizontal incision (1,5-2,0 cm) and blunt preparation up to trachea




Puncture of trachea




Insertion of guide wire (Seldinger wire)



Pre-dilation of injection canal by pre-dilator, removing of pre-dilator



Moisten dilator. Push dilator with inserted guiding catheter into the trachea via guide wire up to the required size of tracheal orifice. After that, remove dilator.



Place the atraumatic inserter in the tracheostomy tube, coat the cuff which has already been checked for leaks with lubricant gel and introduce it into the trachea orifice via the guiding catheter with guide wire inside. After positioning, remove the atraumatic inserter, guiding catheter and guide wire together.

Package contents per set:

PercuTrach® Dilation set I

- Pre-dilator
- Puncture cannula (14G)
- Dilator (hydrophilic coating)
- Guiding catheter with safety lock
- Guide wire with dispenser
- Gauze compresses
- Syringe
- Scalpel

Optima PercuTrach® Basic Cuff


- Outer cannula with low-pressure cuff and 15 mm connector
- Atraumatic inserter
- 1 transparent, unfenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Neck strap
- Lubricant gel



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 1
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112424
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112425
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112426

INCLUDING:



NECK STRAP

OBTURATOR

LUBRICANT GEL

Package contents per set:

PercuTrach® Dilation set I

- Pre-dilator
- Puncture cannula (14G)
- Dilator (hydrophilic coating)
- Guiding catheter with safety lock
- Guide wire with dispenser
- Gauze compresses
- Syringe
- Scalpel

Optima PercuTrach® Voice Cuff


- Sieved outer cannula with low-pres-
sure cuff and 15 mm connector
- Atraumatic inserter
- 1 red, fenestrated inner cannula
- 1 white, unfenestrated, radiopaque
inner cannula
- Obturator, perforated
- Decannulation cap
- Speaking valve
- Neck strap
- Lubricant gel



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 1
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112434
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112435
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112436

INCLUDING:



Package contents per set:

PercuTrach® Dilation set I

- Pre-dilator
- Puncture cannula (14G)
- Dilator (hydrophilic coating)
- Guiding catheter with safety lock
- Guide wire with dispenser
- Gauze compresses
- Syringe
- Scalpel

Optima PercuTrach® Suction


- Outer cannula with low-pressure cuff, subglottic suction function and 15 mm connector
- Atraumatic inserter
- 1 transparent, unfenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Funnel connector
- Neck strap
- Lubricant gel



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 1
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112454
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112455
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112456

INCLUDING:



Package contents per set:

PercuTrach® Dilation set I

- Pre-dilator
- Puncture cannula (14G)
- Dilator (hydrophilic coating)
- Guiding catheter with safety lock
- Guide wire with dispenser
- Gauze compresses
- Syringe
- Scalpel

Optima PercuTrach® Suction Voice

- Sieved outer cannula with low-pressure cuff, subglottic suction function and 15 mm connector
- Atraumatic inserter
- 1 red, fenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Funnel connector
- Decannulation cap
- Speaking valve
- Neck strap
- Lubricant gel



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 1
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112444
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112445
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set I + Optima Percutrach® tracheostomy tube)	112446

INCLUDING:



NECK STRAP



OBTURATOR



LUBRICANT GEL



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



FUNNEL
CONNECTOR

DILATATION TRACHEOSTOMY

developed by Griggs


The percutaneous dilatation tracheostomy technique according to Griggs is a minimal-invasive treatment. This method implies the tracheal dilation by a tracheal dilating forceps. The closed forceps, featuring a longitudinal duct, is pushed into the trachea via a guide wire. Then the forceps is spread and thus dilates the trachea until the insertion of the required tracheostomy tube is possible.

This technique is a quick and efficient, tissue-conserving dilation procedure and, similar to the Ciaglia method, is characterized by a small skin incision and causing few complications.


The procedure should be performed under bronchoscopic control.




Position patient with hyperextension of neck, cover area around the neck with sterile dressing, disinfect skin




Palpation of neck, sonography (recommended), diaphanoscopy of neck area (obligatory), horizontal incision (1,5-2,0 cm) and blunt preparation up to trachea




Puncture of trachea




Insertion of guide wire (Seldinger wire)



Opening of injection canal by pre-dilator




Insert closed tracheal dilating forceps via guide wire and dilate pre-tracheal tissue first.



Insert forceps further into trachea and dilate tracheal tissue by spreading.

Remove spread forceps.



Check the cuff of the tracheostomy tube for leaks. Then introduce the tracheostomy tube with the inserted obturator into the trachea orifice via the guide wire. After positioning, remove the guide wire and obturator together.

Package contents per set:

PercuTrach® Dilation set II

- Tracheal dilating forceps
- Pre-dilator
- Puncture cannula (14G)
- Gauze compresses
- Syringe
- Scalpel
- Guide wire with dispenser

Optima Basic Cuff


- Outer cannula with cuff and 15 mm universal connector
- 1 transparent, unfenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Neck strap



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 2
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112324
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112325
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112326

INCLUDING:



NECK STRAP

OBTURATOR


Package contents per set:

PercuTrach® Dilation set II

- Tracheal dilating forceps
- Pre-dilator
- Puncture cannula (14G)
- Gauze compresses
- Syringe
- Scalpel
- Guide wire with dispenser

Optima Voice Cuff

- Sieved outer cannula with cuff and 15 mm universal connector
- 1 red, fenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Decannulation cap
- Speaking valve
- Neck strap



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 2
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112334
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112335
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112336

INCLUDING:




Package contents per set:

PercuTrach® Dilation set II

- Tracheal dilating forceps
- Pre-dilator
- Puncture cannula (14G)
- Gauze compresses
- Syringe
- Scalpel
- Guide wire with dispenser

Optima Suction


- Outer cannula with cuff, subglottic suction function and 15 mm universal connector
- 1 transparent, unfenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Funnel connector
- Neck strap



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 2
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112354
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112355
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112356

INCLUDING:



NECK STRAP

OBTURATOR

FUNNEL CONNECTOR

Package contents per set:

PercuTrach® Dilation set II

- Tracheal dilating forceps
- Pre-dilator
- Puncture cannula (14G)
- Gauze compresses
- Syringe
- Scalpel
- Guide wire with dispenser

Optima Suction Voice


- Sieved outer cannula with cuff, subglottic suction function and 15 mm universal connector
- 1 red, fenestrated inner cannula
- 1 white, unfenestrated, radiopaque inner cannula
- Obturator, perforated
- Funnel connector
- Decannulation cap
- Speaking valve
- Neck strap



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length along centre line (OC)	Length along outer bend (OC)	PU	Article number PDT SET TYPE 2
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112344
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112345
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 Complete Set (Dilation set II + Optima tracheostomy tube)	112346

INCLUDING:



Sets for percutaneous dilatational
tracheostomy without tracheostomy tube

Dilation set I

- Pre-dilator
- Puncture cannula (14G)
- Dilator (hydrophilic coating)
- Guiding catheter with safety lock
- Guide wire with dispenser
- Gauze compresses
- Syringe
- Scalpel



Dilation set II

- Tracheal dilating forceps
- Pre-dilator
- Puncture cannula (14G)
- Gauze compresses
- Syringe
- Scalpel
- Guide wire with dispenser



	PU	Article number
Dilation set I	1 set	102086
Dilation set II	1 set	102087

Dilation set II can also be used in conjunction with
Optima Percutrach® tracheostomy tubes

Guide wire with dispenser
(Seldinger wire)

- Made of kink-resistant NiTiNOL
- Individually sterile packaged



	PU	Article number
Guide wire with dispenser (Seldinger Wire)	5 pieces	102090



TRACHEOSTOMY TUBES

OPTIMA PercuTrach® tracheostomy tubes with atraumatic inserters have been specially developed for use during a dilatational tracheostomy.

The thermosensitive material of the outer cannula with a rigid 15 mm universal connector ensures pleasant wearing comfort. The atraumatic inserter is perfectly adapted for this purpose and snaps securely into the outer cannula. Moreover, it features a practical handle which incorporates elevations that prevent slipping.


The tracheostomy tube set is always provided with two inner cannulas which can be inserted and removed without difficulty and are fixed securely in the outer cannula.

OPTIMA PercuTrach® tracheostomy tubes always have a thin-walled low-pressure cuff which is filled via a one-way valve and securely seals the cannula against the tracheal wall. This avoids inadvertent inhalation of secretions. The soft cuff adapts superbly to the trachea anatomy and ensures reliable sealing, even where low filling pressure is involved.

OPTIMA PercuTrach® tracheostomy tubes are available in sizes 7 – 9 and different designs.

All materials are biocompatible and free of latex and DEHP.

The 15 mm outer cannula connector enables attachment of accessories (HMEs and speaking valves) and serves as a connection for ventilators,- with or without the use of an inner cannula.



MEASUREMENTS FOR OPTIMA PERCUTRACH®:

Size	ID _{OC}	ID _{IC}	OD	L _n Length along centre line (OC)	L _b Length along outer bend (OC)
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm

ID_{OC} = Inside diameter (clear width) at bottom of outer cannula
ID_{IC} = Inside diameter at bottom of inner cannula
OD = Outside diameter of outer cannula behind the flange
IC = Inner cannula
OC = Outer cannula
L_n = Length along centre line of the outer cannula from neck flange to bottom of the tube
L_b = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit

Product features

- Outer cannula with 15 mm connector and low-pressure cuff
- With atraumatic inserter which ensures gentle and safe positioning of the tracheostomy tube during dilatational tracheostomy or subsequent tube changing
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free



Indication


- Dilatation tracheostomy
- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation

OPTIMA PERCUTRACH®
BASIC CUFF

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 set	102424
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 set	102425
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 set	102426

INCLUDING:



OPTIONAL
ACCESSORIES

Product features

- Sieved outer cannula with 15 mm connector and low-pressure cuff
- With atraumatic inserter which ensures gentle and safe positioning of the tracheostomy tube during dilatational tracheostomy or subsequent tube changing
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free



Indication

- Dilatation tracheostomy
- Tracheostomy
- Phonation
- Ventilation

OPTIMA PERCUTRACH®
VOICE CUFF

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln Lb		PU	Article number
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 set	102434
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 set	102435
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 set	102436

INCLUDING:



NECK STRAP



OBTURATOR



LUBRICANT GEL



PRIMED
SPEAKING VALVE



DECANNULATION
CAP




OPTIONAL
ACCESSORIES

Product features

- Outer cannula with 15 mm connector, low-pressure cuff and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- With atraumatic inserter which ensures gentle and safe positioning of the tracheostomy tube during dilatational tracheostomy or subsequent tube changing
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Dilatation tracheostomy
- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation



OPTIMA PERCUTRACH®
SUCTION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length LnLb		PU	Article number
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 set	102454
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 set	102455
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 set	102456

INCLUDING:



NECK STRAP



OBTURATOR



LUBRICANT GEL



FUNNEL CONNECTOR



OPTIONAL
ACCESSORIES


TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula with 15 mm connector, low-pressure cuff and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- With atraumatic inserter which ensures gentle and safe positioning of the tracheostomy tube during dilatational tracheostomy or subsequent tube changing
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Dilatation tracheostomy
- Tracheostomy
- Phonation
- Ventilation




OPTIMA PERCUTRACH®
SUCTION VOICE


SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln Lb		PU	Article number
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm	1 set	102444
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm	1 set	102445
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm	1 set	102446


INCLUDING:




NECK STRAP




OBTURATOR




LUBRICANT GEL



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



FUNNEL
CONNECTOR



OPTIONAL
ACCESSORIES



TRACHEOSTOMY TUBES

All types of **OPTIMA**® tracheostomy tubes are characterized by a transparent outer cannula with a firmly attached universal 15 mm ISO connector and flexible thermosensitive material properties which ensure good wear properties.


The tubes are lightweight and the inserting and removing of the inner cannula is simple. The inner cannula can be fitted into the 15 mm connector of the outer cannula and is easy to grip via lateral grooves.

The soft flange has lateral eyelets for fixing the inserted tube with a neck strap.

The **OPTIMA**® series comprises a complete range of tracheostomy tubes for tracheostomized patients. They are available in five adult sizes 6-10, as uncuffed or cuffed types, with and without phonation or subglottic suction function.

All materials are biocompatible, latex-free and DEHP-free.

With or without inner cannula in place, the built-in 15 mm ISO connector of the outer cannula allows the attachment of accessories (heat and moisture exchangers (HMEs) and speaking valves) and can also be used for attaching ventilators (in case of cuffed tubes).



MEASUREMENTS FOR OPTIMA:

Size	ID_{OC}	ID_{IC}	OD	Ln Length along centre line (OC)	Lb Length along outer bend (OC)
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73,0 mm
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100,0 mm

ID_{OC} = Inside diameter (clear width) at bottom of outer cannula

ID_{IC} = Inside diameter at bottom of inner cannula

OD = Outside diameter of outer cannula behind the flange

IC = Inner cannula

OC = Outer cannula

Ln = Length along centre line of the outer cannula from neck flange to bottom of the tube

Lb = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit

TRACHEOSTOMY TUBES

Product features

- Outer cannula with 15 mm connector
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI




OPTIMA BASIC

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	102013
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	102014
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	102015
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	102016
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	102017

INCLUDING:



OPTIONAL ACCESSORIES

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula with 15 mm connector
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI



OPTIMA VOICE

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	102033
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	102034
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	102035
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	102036
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	102037

INCLUDING:



NECK STRAP



OBTURATOR



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

Product features

- Outer cannula with 15 mm connector and low-pressure cuff
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation



OPTIMA BASIC CUFF

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	102023
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	102024
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	102025
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	102026
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	102027

INCLUDING:



OPTIONAL ACCESSORIES

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula with 15 mm connector and low-pressure cuff
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy
- Phonation
- Ventilation




OPTIMA VOICE CUFF

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	102043
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	102044
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	102045
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	102046
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	102047


INCLUDING:



NECK STRAP



OBTURATOR



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

Product features

- Outer cannula with 15 mm connector, low-pressure cuff and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation




OPTIMA SUCTION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	102053
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	102054
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	102055
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	102056
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	102057

INCLUDING:




OPTIONAL ACCESSORIES

Product features

- Sieved outer cannula with 15 mm connector, low-pressure cuff and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy
- Phonation
- Ventilation




OPTIMA SUCTION VOICE


SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	102063
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	102064
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	102065
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	102066
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	102067


INCLUDING:




NECK STRAP




OBTURATOR



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



FUNNEL CONNECTOR



OPTIONAL
ACCESSORIES

REPLACEMENT INNER
CANNULAS

Inner cannula A is unfenestrated and transparent for improved secretion assessment

Inner cannula B is red, fenestrated and enables phonation in conjunction with a speaking valve which is connected to the 15 mm connector of the outer cannula

Inner cannula C is unfenestrated, white and radiopaque.

(A)



(B)



(C)





SIZE CHART:

REPLACEMENT INNER CANNULAS				
Size	PU	Transparent (A)	Red, fenestrated (B)	White, radiopaque (C)
		Article number	Article number	Article number
6	5/Box	102018	102028	102038
7	5/Box	102019	102029	102039
8	5/Box	102020	102030	102040
9	5/Box	102021	102031	102041
10	5/Box	102022	102032	102042



ULTRASOFT tracheostomy tubes are distinguished by particularly soft and flexible material. This applies both to the outer cannula and the inner cannulas. This guarantees wearing properties which ensure maximum wearing comfort.

Its multidimensional movable neck flange enables adaptation of the **ULTRASOFT** tracheostomy tube to the individual anatomical conditions of every patient. The head remains freely movable without exerting pressure on the trachea.



Variants with a cuff are equipped with an extremely thin-walled cylindrical low-pressure polyurethane (PU) cuff which adapts ideally to the trachea. This cuff is functionally superior to conventional PVC balloons and minimises microaspiration to the greatest possible extent.


All materials are biocompatible, latex-free and DEHP-free.

The **ULTRASOFT** series comprises a complete range of tracheostomy tubes for tracheostomized patients. They are available in five adult sizes 6-10, as uncuffed or cuffed types, with and without phonation or subglottic suction function.

All types of **ULTRASOFT** tracheostomy tubes are characterized by a white outer cannula with a firmly attached universal 15 mm connector. With or without inner cannula in place, this connector allows the attachment of accessories (heat and moisture exchangers (HME) and speaking valves) and can also be used to attach ventilators (in case of cuffed tubes).

Handling is easy and safe. The inner cannula can be inserted in the outer cannula and fitted into the 15 mm connector. The lateral grooves make it easy to grip and remove the inner cannula.

The lateral eyelets enable the fixing of the inserted tube with a neck strap.



MEASUREMENTS FOR ULTRASOFT:

Size	ID_{OC}	ID_{IC}	OD	Ln Length along centre line (OC)	Lb Length along outer bend (OC)
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73,0 mm
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79,0 mm
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87,0 mm
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93,0 mm
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100,0 mm

ID_{OC} = Inside diameter (clear width) at bottom of outer cannula
 ID_{IC} = Inside diameter at bottom of inner cannula
 OD = Outside diameter of outer cannula behind the flange
 IC = Inner cannula
 OC = Outer cannula
 Ln = Length along centre line of the outer cannula from neck flange to bottom of the tube
 Lb = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit

Product features

- Outer cannula with 15 mm connector
- High wearing comfort due to especially soft material of both outer and inner cannula
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI



ULTRASOFT BASIC

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	103013
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	103014
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	103015
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	103016
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	103017

INCLUDING:



OPTIONAL ACCESSORIES

Product features

- Sieved outer cannula with 15 mm connector
- High wearing comfort due to especially soft material of both outer and inner cannula
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, latex-free, DEHP-free

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI



ULTRASOFT VOICE


SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	103033
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	103034
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	103035
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	103036
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	103037

INCLUDING:



NECK STRAP



OBTURATOR



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

Product features

- Outer cannula with 15 mm connector and particularly thin-walled Polyurethan (PU) low-pressure cuff
- High wearing comfort due to especially soft material of both outer and inner cannula
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, PU, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation




ULTRASOFT BASIC CUFF

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	103023
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	103024
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	103025
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	103026
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	103027

INCLUDING:



OPTIONAL ACCESSORIES

Product features

- Sieved outer cannula with 15 mm connector and particularly thin-walled Polyurethan (PU) low-pressure cuff
- High wearing comfort due to especially soft material of both outer and inner cannula
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, PU, latex-free, DEHP-free

Indication

- Tracheostomy
- Phonation
- Ventilation



ULTRASOFT VOICE CUFF

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	103043
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	103044
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	103045
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	103046
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	103047

INCLUDING:



NECK STRAP



OBTURATOR



PRIMED
SPEAKING VALVE



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

Product features

- Outer cannula with 15 mm connector, particularly thin-walled Polyurethan (PU) low-pressure cuff and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- High wearing comfort due to especially soft material of both outer and inner cannula
- With inner cannula system
- 1 transparent inner cannula for improved secretion assessment
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, PU, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation




ULTRASOFT SUCTION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	103053
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	103054
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	103055
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	103056
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	103057

INCLUDING:



OPTIONAL ACCESSORIES

Product features

- Sieved outer cannula with 15 mm connector, particularly thin-walled Polyurethan (PU) low-pressure cuff and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- High wearing comfort due to especially soft material of both outer and inner cannula
- With inner cannula system
- 1 red, fenestrated inner cannula enables phonation in conjunction with a speaking valve
- 1 white inner cannula made of radiopaque material for position control
- Material: PVC, PU, latex-free, DEHP-free

Indication

- Tracheostomy
- Phonation
- Ventilation




ULTRASOFT SUCTION VOICE

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,0 mm	4,4 mm	9,2 mm	64,5 mm	73 mm	1 set	103063
7	7,0 mm	5,4 mm	10,5 mm	70,0 mm	79 mm	1 set	103064
8	8,0 mm	6,4 mm	11,6 mm	75,5 mm	87 mm	1 set	103065
9	9,0 mm	7,4 mm	12,6 mm	81,0 mm	93 mm	1 set	103066
10	10,0 mm	8,4 mm	13,6 mm	87,5 mm	100 mm	1 set	103067

INCLUDING:



OPTIONAL ACCESSORIES

REPLACEMENT INNER
CANNULAS

Inner cannula A is unfenestrated and transparent for improved secretion assessment

(A)



Inner cannula B is red, fenestrated and enables phonation in conjunction with a speaking valve which is connected to the 15 mm connector of the outer cannula

(B)



Inner cannula C is unfenestrated, white and radiopaque.

(C)



SIZE CHART:

REPLACEMENT INNER CANNULAS				
		Transparent (A)	Red, fenestrated (B)	White, radiopaque (C)
Size	PU	Article number	Article number	Article number
6	5/Box	102018	102028	102038
7	5/Box	102019	102029	102039
8	5/Box	102020	102030	102040
9	5/Box	102021	102031	102041
10	5/Box	102022	102032	102042



TRACHEOSTOMY TUBES

PRIFLEX® tracheostomy tubes are spiral-reinforced cannulas that combine a highly stable shape with soft and flexible material properties.

The flexible **PRIFLEX®** tracheostomy tubes feature an individually adjustable neck flange with soft and flexible wings which ensures accurate placement of the cannula. A scale on the tracheostomy tube facilitates the exact positioning.

PRIFLEX® tracheostomy tubes feature a firm 15 mm universal connector which enables the attachment of accessories (HME and speaking valves) as well as mechanical ventilators.


The metal spiral in the wall of the soft tube ensures that the tube retains its shape and it prevents kinking. Even the **PRIFLEX®** tracheostomy tubes with speech function have a continuous spiral reinforcement, in which a phonation sieving is imbedded.



Sieving in continuous reinforcement
(PRIFLEX® tracheostomy tube)

PRIFLEX® tracheostomy tubes are high-quality sterile products made of latex-free medical-grade materials.

These tracheostomy tubes are appropriate for insertion especially in ventilated patients, adipose patients and in case of postoperative oedema as well as a sunken tracheostoma.



MEASUREMENTS FOR PRIFLEX®:

for the types

PRIFLEX® without Cuff, PRIFLEX® with Cuff and PRIFLEX® Suction


Size	ID	OD	L Lenght along centre line adjustable to
6	6,0 mm	9,4 mm	69,0 mm
7	7,0 mm	10,4 mm	69,0 mm
8	8,0 mm	11,4 mm	96,0 mm
9	9,0 mm	12,4 mm	122,0 mm
10	10,0 mm	13,4 mm	122,0 mm
11	11,0 mm	14,4 mm	123,0 mm

ID = Inside diameter at the bottom of the tube

OD = Outside diameter at the bottom of the tube

L = Length along centre line of the tube from neck flange to bottom of the tube

PU = Packaging unit



MEASUREMENTS FOR PRIFLEX® XL:

Size	ID	OD	L Lenght along centre line adjustable to
7	7,0 mm	10,4 mm	155,0 mm
8	8,0 mm	11,4 mm	155,0 mm
9	9,0 mm	12,4 mm	155,0 mm
10	10,0 mm	13,4 mm	155,0 mm
11	11,0 mm	14,4 mm	155,0 mm

ID = Inside diameter at the bottom of the tube

OD = Outside diameter at the bottom of the tube

L = Length along centre line of the tube from neck flange to bottom of the tube

PU = Packaging unit

TRACHEOSTOMY TUBES

Product features

- Spiral-reinforced tracheostomy tube with 15 mm connector and individually adjustable neck flange
- Versions with speaking option feature a special kink-resistant sieving within a continuous spiral-reinforcement
- Material: PVC, metal spiral, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Obesity/ Oedema



WITHOUT CUFF



WITHOUT CUFF/ WITH SIEVING


SIZE CHART:

Size	ID	OD	Length (adjustable to)	PU	WITHOUT CUFF	WITHOUT CUFF/ WITH SIEVING
					Article number	Article number
6	6,0 mm	9,4 mm	69 mm	1 set	200802	200803
7	7,0 mm	10,4 mm	69 mm	1 set	200804	200809
8	8,0 mm	11,4 mm	96 mm	1 set	200805	200810
9	9,0 mm	12,4 mm	122 mm	1 set	200806	200811
10	10,0 mm	13,4 mm	122 mm	1 set	200807	200812
11	11,0 mm	14,4 mm	123 mm	1 set	200808	200813

INCLUDING:



NECK STRAP



OBTURATOR
WITH SUCTION FUNCTION



OPTIONAL
ACCESSORIES

Product features

- Spiral-reinforced tracheostomy tube with 15 mm connector, low-pressure cuff and individually adjustable neck flange
- Versions with speaking option feature a special kink-resistant sieving within a continuous spiral-reinforcement
- Material: PVC, metal spiral, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Ventilation
- Obesity/ Oedema



WITH CUFF



WITH CUFF/ WITH SIEVING


SIZE CHART:

Size	ID	OD	Length (adjustable to)	PU	WITHOUT CUFF	WITHOUT CUFF/ WITH SIEVING
					Article number	Article number
6	6,0 mm	9,4 mm	69 mm	1 set	200824	200825
7	7,0 mm	10,4 mm	69 mm	1 set	200814	200819
8	8,0 mm	11,4 mm	96 mm	1 set	200815	200820
9	9,0 mm	12,4 mm	122 mm	1 set	200816	200821
10	10,0 mm	13,4 mm	122 mm	1 set	200817	200822
11	11,0 mm	14,4 mm	123 mm	1 set	200818	200823

INCLUDING:



NECK STRAP



OBTURATOR
WITH SUCTION FUNCTION



OPTIONAL
ACCESSORIES

Product features

- Spiral-reinforced tracheostomy tube with 15 mm connector, low-pressure cuff, individually adjustable neck flange and subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand)
- Material: PVC, metal spiral, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation
- Obesity/ Oedema



PRIFLEX® SUCTION

SIZE CHART:

Size	ID	OD	Length (adjustable to)	PU	Article number
6	6,0 mm	9,4 mm	69 mm	1 set	203713
7	7,0 mm	10,4 mm	69 mm	1 set	203714
8	8,0 mm	11,4 mm	96 mm	1 set	203715
9	9,0 mm	12,4 mm	122 mm	1 set	203716
10	10,0 mm	13,4 mm	122 mm	1 set	203717
11	11,0 mm	14,4 mm	123 mm	1 set	203718

INCLUDING:



NECK STRAP



OBTURATOR
WITH SUCTION FUNCTION



FUNNEL CONNECTOR



OPTIONAL
ACCESSORIES

Product features

- Extra-long spiral-reinforced tracheostomy tube with 15 mm connector, low-pressure cuff and individually adjustable neck flange
- Not performed, thus enables individual adaptation to the anatomical conditions of the patient
- Material: PVC, metal spiral, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation
- Obesity/ Oedema
- Tracheal stenoses
- Presence of large parastomal soft tissue (e.g. tumours)



PRIFLEX® XL


SIZE CHART:

Size	ID	OD	Length (adjustable to)	PU	Article number
7	7,0 mm	10,4 mm	155 mm	1 set	200797
8	8,0 mm	11,4 mm	155 mm	1 set	200798
9	9,0 mm	12,4 mm	155 mm	1 set	200799
10	10,0 mm	13,4 mm	155 mm	1 set	200800
11	11,0 mm	14,4 mm	155 mm	1 set	200801

INCLUDING:



NECK STRAP



OBTURATOR
WITH SUCTION FUNCTION



OPTIONAL
ACCESSORIES



TRACHEOSTOMY TUBES

PRIFLEX® MRI stand for a new generation of tracheostomy tubes with plastic spiral reinforcement which are suitable for MRI.


These tracheostomy tubes **do not contain any metal parts** - even the pilot balloon is equipped with a special plastic spring -, what makes them suitable for patients who undergo a magnetic resonance imaging (MRI) procedure or need radiation therapy.

They combine a high form stability with soft and flexible material features, which is also characteristic for tubes with classical metal spiral reinforcement. The soft and flexible tube flange is variably adjustable by a simple screw cap. It adapts exactly to the individual patient requirements. A scale on the tracheostomy tube facilitates the exact positioning, and a firm 15 mm universal connector enables the attachment of accessories (HME and speaking valves) as well as mechanical ventilators.

The integrated plastic spiral ensures the form stability and kink resistance of the tracheostomy tubes.

PRIFLEX® MRI tracheostomy tubes are suitable especially for the use in ventilated neurological patients, as well as in patients whose treatment plan includes MRI or radiation therapy.

TRACHEOSTOMY TUBES



MEASUREMENTS FOR PRIFLEX® MRI:

Size	ID	OD	L Length along centre line adjustable to
6	6,0 mm	8,9 mm	70,0 mm
7	7,0 mm	9,9 mm	80,0 mm
8	8,0 mm	10,9 mm	100,0 mm
9	9,0 mm	11,9 mm	120,0 mm
10	10,0 mm	12,9 mm	130,0 mm
11	11,0 mm	13,9 mm	130,0 mm

ID = Inside diameter at the bottom of the tube

OD = Outside diameter at the bottom of the tube

L = Length along centre line of the tube from neck flange to bottom of the tube

PU = Packaging unit

TRACHEOSTOMY TUBES

Product features

- Metal-free tracheostomy tube with plastic spiral reinforcement, with 15 mm connector and individually adjustable neck flange
- Material: PVC, plastic spiral, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Obesity/ Oedema
- Radiation/ MRI



SIZE CHART:

Size	ID	OD	Length (adjustable to)	PU	Article number
6	6,0 mm	8,9 mm	70 mm	1 set	203802
7	7,0 mm	9,9 mm	80 mm	1 set	203804
8	8,0 mm	10,9 mm	100 mm	1 set	203805
9	9,0 mm	11,9 mm	120 mm	1 set	203806
10	10,0 mm	12,9 mm	130 mm	1 set	203807
11	11,0 mm	13,9 mm	130 mm	1 set	203808

INCLUDING:



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES

Product features

- Metal-free tracheostomy tube with plastic spiral reinforcement, with 15 mm connector, low-pressure cuff and individually adjustable neck flange
- Material: PVC, plastic spiral, latex-free, DEHP-free

Indication

- Tracheostomy/ Laryngectomy
- Ventilation
- Obesity/ Oedema
- Radiation/ MRI



SIZE CHART:

Size	ID	OD	Length (adjustable to)	PU	Article number
6	6,0 mm	8,9 mm	70 mm	1 set	203824
7	7,0 mm	9,9 mm	80 mm	1 set	203814
8	8,0 mm	10,9 mm	100 mm	1 set	203815
9	9,0 mm	11,9 mm	120 mm	1 set	203816
10	10,0 mm	12,9 mm	130 mm	1 set	203817
11	11,0 mm	13,9 mm	130 mm	1 set	203818

INCLUDING:



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS



PRIMEDISTOM® cuffless tracheostomy tubes are designed for permanent wearers after tracheostomy and laryngectomy.


Their soft, flexible and extremely lightweight material guarantees maximum wearing comfort. **PRIMEDISTOM®** tracheostomy tubes are thermo-sensitive and reach their ideal properties at body temperature. Due to these material characteristics **PRIMEDISTOM®** tracheostomy tubes adapt very well to the shape of the trachea.

All materials are biocompatible, latex-free and DEHP-free.

PRIMEDISTOM® cuffless tracheostomy tubes are offered in a wide range of different types and sizes to ensure that all individual patient's requirements are met.

Customized service is available on request.

TRACHEOSTOMY TUBES
CUFFLESS



MEASUREMENTS FOR PRIMEDISTOM® CUFFLESS:

Size	ID _{OC}	ID _{IC}	OD	Ln Length along centre line (OC)			Lb Length along outer bend (OC)		
				short version	standard version	long version	short version	standard version	long version
3,5	3,6 mm	2,7 mm	6,2 mm		42 mm			46 mm	
4	4,4 mm	3,3 mm	7,2 mm		50 mm			55 mm	
5	5,4 mm	4,3 mm	8,2 mm		52 mm	87 mm		57 mm	92 mm
6	6,2 mm	4,9 mm	9,2 mm		54 mm	88 mm		60 mm	95 mm
7	7,0 mm	5,8 mm	10,2 mm	40 mm	57 mm	93 mm	44 mm	65 mm	100 mm
8	8,0 mm	6,8 mm	11,2 mm	41 mm	61 mm	97 mm	44 mm	70 mm	105 mm
9	9,0 mm	7,8 mm	12,2 mm	41 mm	65 mm	101 mm	44 mm	75 mm	110 mm
10	10,0 mm	8,8 mm	13,2 mm	43 mm	74 mm	105 mm	44 mm	85 mm	115 mm
11	11,0 mm	9,8 mm	14,2 mm	42 mm	76 mm	109 mm	44 mm	87 mm	120 mm
12	12,0 mm	10,8 mm	15,2 mm	42 mm	78 mm	113 mm	44 mm	90 mm	125 mm
13	13,0 mm	11,2 mm	16,2 mm	42 mm	88 mm	118 mm	44 mm	100 mm	130 mm

ID_{OC} = Inside diameter (clear width) at bottom of outer cannula
ID_{IC} = Inside diameter at bottom of inner cannula
OD = Outside diameter of outer cannula behind the flange
IC = Inner cannula
OC = Outer cannula
Ln = Length along centre line of the outer cannula from neck flange to bottom of the tube
Lb = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Tracheostomy tube with integrated X-ray contrast line for position control
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI



SHORT VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number
			Ln	Lb		
3,5	3,6 mm	6,2 mm	42 mm	46 mm	1 set	200010
4	4,4 mm	7,2 mm	50 mm	55 mm	1 set	200011
5	5,4 mm	8,2 mm	52 mm	57 mm	1 set	201012
6	6,2 mm	9,2 mm	54 mm	60 mm	1 set	201013
7	7,0 mm	10,2 mm	57 mm	65 mm	1 set	201014
8	8,0 mm	11,2 mm	61 mm	70 mm	1 set	201015
9	9,0 mm	12,2 mm	65 mm	75 mm	1 set	201016
10	10,0 mm	13,2 mm	74 mm	85 mm	1 set	201017
11	11,0 mm	14,2 mm	76 mm	87 mm	1 set	201018
12	12,0 mm	15,2 mm	78 mm	90 mm	1 set	201019
13	13,0 mm	16,2 mm	88 mm	100 mm	1 set	201020

SHORT VERSION			
Length		PU	Article number
Ln	Lb		
40 mm	44 mm	1 set	200204
41 mm	44 mm	1 set	200205
41 mm	44 mm	1 set	200206
43 mm	44 mm	1 set	200207
42 mm	44 mm	1 set	200208
42 mm	44 mm	1 set	200209
42 mm	44 mm	1 set	200210

INCLUDING:



NECK STRAP

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Tracheostomy tube with 15 mm connector and integrated X-ray contrast line for position control
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number
			Ln	Lb		
5	5,4 mm	8,2 mm	52 mm	57 mm	1 set	201012
6	6,2 mm	9,2 mm	54 mm	60 mm	1 set	201013
7	7,0 mm	10,2 mm	57 mm	65 mm	1 set	201014
8	8,0 mm	11,2 mm	61 mm	70 mm	1 set	201015
9	9,0 mm	12,2 mm	65 mm	75 mm	1 set	201016
10	10,0 mm	13,2 mm	74 mm	85 mm	1 set	201017
11	11,0 mm	14,2 mm	76 mm	87 mm	1 set	201018
12	12,0 mm	15,2 mm	78 mm	90 mm	1 set	201019
13	13,0 mm	16,2 mm	88 mm	100 mm	1 set	201020

LONG VERSION				
Length	Ln	Lb	PU	Article number
87 mm	92 mm	1 set		201023
88 mm	95 mm	1 set		201024
93 mm	100 mm	1 set		201025
97 mm	105 mm	1 set		201026
101 mm	110 mm	1 set		201027
105 mm	115 mm	1 set		201028
109 mm	120 mm	1 set		201029
113 mm	125 mm	1 set		201030
118 mm	130 mm	1 set		201031

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Tracheostomy tube with 22 mm connector and integrated X-ray contrast line for position control
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number
			Ln	Lb		
5	5,4 mm	8,2 mm	52 mm	57 mm	1 set	201034
6	6,2 mm	9,2 mm	54 mm	60 mm	1 set	201035
7	7,0 mm	10,2 mm	57 mm	65 mm	1 set	201036
8	8,0 mm	11,2 mm	61 mm	70 mm	1 set	201037
9	9,0 mm	12,2 mm	65 mm	75 mm	1 set	201038
10	10,0 mm	13,2 mm	74 mm	85 mm	1 set	201039
11	11,0 mm	14,2 mm	76 mm	87 mm	1 set	201040
12	12,0 mm	15,2 mm	78 mm	90 mm	1 set	201041
13	13,0 mm	16,2 mm	88 mm	100 mm	1 set	201042

LONG VERSION			
Length	Ln	Lb	Article number
PU	87 mm	92 mm	201045
	88 mm	95 mm	201046
1 set	93 mm	100 mm	201047
	97 mm	105 mm	201048
1 set	101 mm	110 mm	201049
	105 mm	115 mm	201050
1 set	109 mm	120 mm	201051
	113 mm	125 mm	201052
1 set	118 mm	130 mm	201053

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved tracheostomy tube with 22 mm connector and integrated X-ray contrast line for position control
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number
			Ln	Lb		
5	5,4 mm	8,2 mm	52 mm	57 mm	1 set	201076
6	6,2 mm	9,2 mm	54 mm	60 mm	1 set	201077
7	7,0 mm	10,2 mm	57 mm	65 mm	1 set	201078
8	8,0 mm	11,2 mm	61 mm	70 mm	1 set	201079
9	9,0 mm	12,2 mm	65 mm	75 mm	1 set	201080
10	10,0 mm	13,2 mm	74 mm	85 mm	1 set	201081
11	11,0 mm	14,2 mm	76 mm	87 mm	1 set	201082
12	12,0 mm	15,2 mm	78 mm	90 mm	1 set	201083
13	13,0 mm	16,2 mm	88 mm	100 mm	1 set	201084

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	201056
88 mm	95 mm	1 set	201057
93 mm	100 mm	1 set	201058
97 mm	105 mm	1 set	201059
101 mm	110 mm	1 set	201060
105 mm	115 mm	1 set	201061
109 mm	120 mm	1 set	201062
113 mm	125 mm	1 set	201063
118 mm	130 mm	1 set	201064

INCLUDING:



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Short version of tracheostomy tube see p. 68



SIZE CHART:

Size	ID	SHORT VERSION			PU	Article number
		OD	Length			
			Ln	Lb		
9	9,0 mm	12,2 mm	41 mm	44 mm	1 set	200500
10	10,0 mm	13,2 mm	43 mm	44 mm	1 set	200501
11	11,0 mm	14,2 mm	42 mm	44 mm	1 set	200502
12	12,0 mm	15,2 mm	42 mm	44 mm	1 set	200503
13	13,0 mm	16,2 mm	42 mm	44 mm	1 set	200504

INCLUDING:



TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Outer cannula with integrated X-ray contrast line for position control
- 1 inner cannula without connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
3,5	3,6 mm	2,7 mm	6,2 mm	42 mm	46 mm	1 set	200110
4	4,4 mm	3,3 mm	7,2 mm	50 mm	55 mm	1 set	200111
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	200112
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	200113
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	200114
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	200115
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	200116
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	200117
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200118
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	200119
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	200120

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	200122
88 mm	95 mm	1 set	200123
93 mm	100 mm	1 set	200124
97 mm	105 mm	1 set	200125
101 mm	110 mm	1 set	200126
105 mm	115 mm	1 set	200127
109 mm	120 mm	1 set	200128
113 mm	125 mm	1 set	200129
118 mm	130 mm	1 set	200130

INCLUDING:



NECK STRAP



EASY-CLICK

A fixing ring is attached to the inner cannula which is connected through engaging with the outer cannula. Furthermore, the inner cannula can be removed again with ease (pull forwards-downwards without rotation!)

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Outer cannula with integrated X-ray contrast line for position control
- 1 inner cannula with 15 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART: :

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201112
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201113
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201114
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201115
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201116
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201117
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201118
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201119
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201120

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	201123
88 mm	95 mm	1 set	201124
93 mm	100 mm	1 set	201125
97 mm	105 mm	1 set	201126
101 mm	110 mm	1 set	201127
105 mm	115 mm	1 set	201128
109 mm	120 mm	1 set	201129
113 mm	125 mm	1 set	201130
118 mm	130 mm	1 set	201131

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula without connector
- 1 inner cannula with 15 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
3,5	3,6 mm	2,7 mm	6,2 mm	42 mm	46 mm	1 set	200400
4	4,4 mm	3,3 mm	7,2 mm	50 mm	55 mm	1 set	200401
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	200402
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	200403
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	200404
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	200405
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	200406
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	200407
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200408
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	200409
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	200410

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	200412
88 mm	95 mm	1 set	200413
93 mm	100 mm	1 set	200414
97 mm	105 mm	1 set	200415
101 mm	110 mm	1 set	200416
105 mm	115 mm	1 set	200417
109 mm	120 mm	1 set	200418
113 mm	125 mm	1 set	200419
118 mm	130 mm	1 set	200420

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 2 inner cannulas with 15 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201242
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201243
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201244
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201245
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201246
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201247
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201248
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201249
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201250

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula with 15 mm connector
- 1 inner cannula with 22 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201212
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201213
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201214
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201215
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201216
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201217
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201218
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201219
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201220

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	201273
88 mm	95 mm	1 set	201274
93 mm	100 mm	1 set	201275
97 mm	105 mm	1 set	201276
101 mm	110 mm	1 set	201277
105 mm	115 mm	1 set	201278
109 mm	120 mm	1 set	201279
113 mm	125 mm	1 set	201280
118 mm	130 mm	1 set	201281

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula without connector
- 1 inner cannula with 22 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:


Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201262
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201263
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201264
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201265
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201266
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201267
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201268
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201269
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201270

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	201223
88 mm	95 mm	1 set	201224
93 mm	100 mm	1 set	201225
97 mm	105 mm	1 set	201226
101 mm	110 mm	1 set	201227
105 mm	115 mm	1 set	201228
109 mm	120 mm	1 set	201229
113 mm	125 mm	1 set	201230
118 mm	130 mm	1 set	201231

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula without connector
- 1 fenestrated inner cannula with silver speaking valve enables phonation
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, sterling silver (925), latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201302
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201303
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201304
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201305
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201306
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201307
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201308
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201309
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201310

LONG VERSION			
Length	Ln	Lb	Article number
87 mm	92 mm	1 set	201383
88 mm	95 mm	1 set	201384
93 mm	100 mm	1 set	201385
97 mm	105 mm	1 set	201386
101 mm	110 mm	1 set	201387
105 mm	115 mm	1 set	201388
109 mm	120 mm	1 set	201389
113 mm	125 mm	1 set	201390
118 mm	130 mm	1 set	201391

INCLUDING:



NECK STRAP



DECANNULATION
CAP

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 fenestrated inner cannula with silver speaking valve enables phonation
- 1 inner cannula with 15 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, sterling silver (925), latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
4	4,4 mm	3,3 mm	7,2 mm	50 mm	55 mm	1 set	200301
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	200302
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	200303
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	200304
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	200305
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	200306
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	200307
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200308
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	200309
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	200310

LONG VERSION			
Length Ln	Length Lb	PU	Article number
87 mm	92 mm	1 set	200322
88 mm	95 mm	1 set	200323
93 mm	100 mm	1 set	200324
97 mm	105 mm	1 set	200325
101 mm	110 mm	1 set	200326
105 mm	115 mm	1 set	200327
109 mm	120 mm	1 set	200328
113 mm	125 mm	1 set	200329
118 mm	130 mm	1 set	200330

INCLUDING:



NECK STRAP



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 fenestrated inner cannula without connector (compatible with Silicone speaking valve to enable phonation)
- 1 inner cannula with 15 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	200343
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	200344
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	200345
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	200346
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	200347
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200348
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	200349
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	200341

INCLUDING:



NECK STRAP



DECANNULATION CAP



SILICONE SPEAKING VALVE



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 fenestrated inner cannula with 15 mm connector enables phonation in conjunction with a speaking valve
- 1 inner cannula with 15 mm connector
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201372
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201373
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201374
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201375
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201376
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201377
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201378
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201379
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201380

INCLUDING:



TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula without connector
- 1 fenestrated inner cannula with 22 mm connector enables phonation in conjunction with a speaking valve
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number	LONG VERSION			
				Ln	Lb			Length Ln	Lb	PU	Article number
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201312	87 mm	92 mm	1 set	201323
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201313	88 mm	95 mm	1 set	201324
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201314	93 mm	100 mm	1 set	201325
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201315	97 mm	105 mm	1 set	201326
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201316	101 mm	110 mm	1 set	201327
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201317	105 mm	115 mm	1 set	201328
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201318	109 mm	120 mm	1 set	201329
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201319	113 mm	125 mm	1 set	201330
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201320	118 mm	130 mm	1 set	201331

INCLUDING:



NECK STRAP



DECANNULATION CAP



PRIMEDI-PHON



TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula without connector
- 1 inner cannula with 15 mm connector
- 1 fenestrated inner cannula with 22 mm connector enables phonation in conjunction with a speaking valve
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Radiation/ MRI
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201342
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201343
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201344
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201345
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201346
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201347
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201348
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201349
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201350

LONG VERSION			
Length	Ln	Lb	PU
87 mm	92 mm	1 set	201353
88 mm	95 mm	1 set	201354
93 mm	100 mm	1 set	201355
97 mm	105 mm	1 set	201356
101 mm	110 mm	1 set	201357
105 mm	115 mm	1 set	201358
109 mm	120 mm	1 set	201359
113 mm	125 mm	1 set	201360
118 mm	130 mm	1 set	201361

INCLUDING:



NECK STRAP



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
CUFFLESS

Product features

- Sieved outer cannula with integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula without connector
- 1 inner cannula with 15 mm connector
- 1 fenestrated inner cannula with silver speaking valve enables phonation
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, sterling silver (925), latex-free, DEHP-free, thermo-sensitive

Indication

- Tracheostomy
- Phonation
- Tracheal stenoses/ Granuloma (for long version)



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
5	5,4 mm	4,3 mm	8,2 mm	52 mm	57 mm	1 set	201535
6	6,2 mm	4,9 mm	9,2 mm	54 mm	60 mm	1 set	201536
7	7,0 mm	5,8 mm	10,2 mm	57 mm	65 mm	1 set	201537
8	8,0 mm	6,8 mm	11,2 mm	61 mm	70 mm	1 set	201538
9	9,0 mm	7,8 mm	12,2 mm	65 mm	75 mm	1 set	201539
10	10,0 mm	8,8 mm	13,2 mm	74 mm	85 mm	1 set	201540
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	201541
12	12,0 mm	10,8 mm	15,2 mm	78 mm	90 mm	1 set	201542
13	13,0 mm	11,2 mm	16,2 mm	88 mm	100 mm	1 set	201543

LONG VERSION			
Length	Ln	Lb	PU
87 mm	92 mm	1 set	201545
88 mm	95 mm	1 set	201546
93 mm	100 mm	1 set	201547
97 mm	105 mm	1 set	201548
101 mm	110 mm	1 set	201549
105 mm	115 mm	1 set	201550
109 mm	120 mm	1 set	201551
113 mm	125 mm	1 set	201552
118 mm	130 mm	1 set	201553

INCLUDING:



NECK STRAP



DECANNULATION
CAP



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
WITH CUFF



TRACHEOSTOMY TUBES WITH CUFF

PRIMEDISTOM® cuffed tracheostomy tubes are made of medical-grade polymer. They are pleasantly flexible and ensure high wearing comfort. The material used is thermosensitive and reaches its ideal properties at body temperature.

The high volume low-pressure cuff, with thin walls and large volume, can be filled through a one-way valve and securely seals the cannula to the tracheal wall. This prevents an accidental aspiration of secretions. The soft cuff fits perfectly into the trachea and ensures a reliable seal, even with a low inflation pressure. The cuff's pilot line is integrated in the outer cannula, achieving a rounded and smooth outer shape. This prevents additional irritation of the mucous membranes.

The key features of this cannula are its easy handling and comfortable use.

PRIMEDISTOM® cuffed tracheostomy tubes are single-patient products, available in different versions and supplied in sterile packaging.


The cuffed **PRIMEDISTOM®** tracheostomy tubes with sieving are used after a tracheostomy with a complete or partially retained larynx and enable the user to speak. They are available with a silicone speaking valve which can simply be attached to the fenestrated inner cannula without connector. In the valve, a flap opens on inspiration allowing air to flow in. The exhaled airflow causes the valve flap to close, thereby enabling phonation.

The cuffed **PRIMEDISTOM®** tracheostomy tubes with subglottic suction function features an integrated suction line for secretion suctioning above the cuff.

All materials are biocompatible, latex-free and DEHP-free.

In addition to the standard versions, **PRIMEDISTOM®** cuffed tracheostomy tubes are also available as custom-made product versions after consultation with the treating physician. In these adapted cannulas, individual anatomical features can be considered (as far as technical possible).

TRACHEOSTOMY TUBES
WITH CUFF



MEASUREMENTS FOR PRIMEDISTOM® WITH CUFF:

Size	ID _{OC}	ID _{IC}	OD	Ln Length along centre line (OC)		Lb Length along outer bend (OC)	
				standard version	long version	standard version	long version
8	8,0 mm	6,8 mm	11,2 mm	62 mm	97 mm	70 mm	105 mm
9	9,0 mm	7,8 mm	12,2 mm	66 mm	96 mm	75 mm	105 mm
10	10,0 mm	8,8 mm	13,2 mm	76 mm	100 mm	85 mm	110 mm
11	11,0 mm	9,8 mm	14,2 mm	76 mm	104 mm	87 mm	115 mm

ID_{OC} = Inside diameter (clear width) at bottom of outer cannula

ID_{IC} = Inside diameter at bottom of inner cannula

OD = Outside diameter of outer cannula behind the flange

IC = Inner cannula

OC = Outer cannula

L_n = Length along centre line of the outer cannula from neck flange to bottom of the tube

L_b = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit

TRACHEOSTOMY TUBES
WITH CUFF

Product features

- Outer cannula with low-pressure cuff and integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula with 15 mm connector
- 1 inner cannula with swivelling 15 mm connector (during ventilation inner cannula has to be locked with the Clip Adapter to prevent it from shifting or removing)
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitiv

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation
- Tracheal stenoses/ Granuloma (for long version)



PRIMEDICUFF®




LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number	LONG VERSION			
				Ln	Lb			Length Ln	Lb	PU	Article number
8	8,0 mm	6,8 mm	11,2 mm	62 mm	70 mm	1 set	200700	97 mm	105 mm	1 set	200704
9	9,0 mm	7,8 mm	12,2 mm	66 mm	75 mm	1 set	200701	96 mm	105 mm	1 set	200705
10	10,0 mm	8,8 mm	13,2 mm	76 mm	85 mm	1 set	200702	100 mm	110 mm	1 set	200706
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200703	104 mm	115 mm	1 set	200707

INCLUDING:



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
WITH CUFF

Product features

- Sieved outer cannula with low-pressure cuff and integrated X-ray contrast line for position control
- With inner cannula system
- 1 fenestrated inner cannula without connector (compatible with Silicone speaking valve)
- 1 inner cannula with swivelling 15 mm connector (during ventilation inner cannula has to be locked with the ClipAdapter to prevent it from shifting or removing)
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitiv

Indication

- Tracheostomy
- Phonation
- Ventilation
- Tracheal stenoses/ Granuloma (for long version)



PRIMEDIPHON®



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
8	8,0 mm	6,8 mm	11,2 mm	62 mm	70 mm	1 set	200710
9	9,0 mm	7,8 mm	12,2 mm	66 mm	75 mm	1 set	200711
10	10,0 mm	8,8 mm	13,2 mm	76 mm	85 mm	1 set	200712
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200713

LONG VERSION			
Length	Ln	Lb	PU Article number
97 mm	105 mm	1 set	200714
96 mm	105 mm	1 set	200715
100 mm	110 mm	1 set	200716
104 mm	115 mm	1 set	200717

INCLUDING:



NECK STRAP



CLIP ADAPTER



SILICONE
SPEAKING VALVE



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
WITH CUFF

Product features

- Outer cannula with low-pressure cuff, subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand) and integrated X-ray contrast line for position control
- With inner cannula system
- 1 inner cannula with 15 mm connector
- 1 inner cannula with swivelling 15 mm connector (during ventilation inner cannula has to be locked with the Clip Adapter to prevent it from shifting or removing)
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitiv

Indication

- Tracheostomy/ Laryngectomy
- Aspiration
- Ventilation
- Tracheal stenoses/ Granuloma (for long version)

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
8	8,0 mm	6,8 mm	11,2 mm	62 mm	70 mm	1 set	200720
9	9,0 mm	7,8 mm	12,2 mm	66 mm	75 mm	1 set	200721
10	10,0 mm	8,8 mm	13,2 mm	76 mm	85 mm	1 set	200722
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200723



PRIMEDISUCTION®



LONG VERSION

LONG VERSION			
Length		PU	Article number
Ln	Lb		
97 mm	105 mm	1 set	200724
96 mm	105 mm	1 set	200725
100 mm	110 mm	1 set	200726
104 mm	115 mm	1 set	200727

INCLUDING:



NECK STRAP



CLIP ADAPTER



FUNNEL CONNECTOR



OPTIONAL
ACCESSORIES

TRACHEOSTOMY TUBES
WITH CUFF

Product features

- Sieved outer cannula with low-pressure cuff, subglottic suction function (features an integrated suction line for subglottic secretion suctioning on demand) and integrated X-ray contrast line for position control
- With inner cannula system
- 1 fenestrated inner cannula without connector (compatible with Silicone speaking valve)
- 1 inner cannula with swivelling 15 mm connector (during ventilation inner cannula has to be locked with the Clip Adapter to prevent it from shifting or removing)
- Secure fixing of inner cannula to outer cannula with Easy-Click-System (see p. 70)
- Conical form facilitates insertion
- Material: PVC, latex-free, DEHP-free, thermo-sensitiv

Indication

- Tracheostomy
- Phonation
- Ventilation
- Tracheal stenoses/ Granuloma (for long version)



PRIMEDICOMB®



LONG VERSION

SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
8	8,0 mm	6,8 mm	11,2 mm	62 mm	70 mm	1 set	200730
9	9,0 mm	7,8 mm	12,2 mm	66 mm	75 mm	1 set	200731
10	10,0 mm	8,8 mm	13,2 mm	76 mm	85 mm	1 set	200732
11	11,0 mm	9,8 mm	14,2 mm	76 mm	87 mm	1 set	200733

LONG VERSION			
Length Ln	Length Lb	PU	Article number
97 mm	105 mm	1 set	200734
96 mm	105 mm	1 set	200735
100 mm	110 mm	1 set	200736
104 mm	115 mm	1 set	200737

INCLUDING:



NECK STRAP



CLIP ADAPTER



SILICONE SPEAKING VALVE



FUNNEL CONNECTOR



OPTIONAL
ACCESSORIES

REPLACEMENT INNER
CANNULAS

Unfenestrated replacement inner cannulas
for PRIMEDISTOM® tracheostomy tube
in standard length

(A)



(B)



(C)



SIZE CHART:

Size	Suitable for Primedistom® outer cannula with Length Lb	PU	Inner cannula unfenestrated, without connector (A) Article number	Inner cannula unfenestrated, with 15 mm connector (B) Article number	Inner cannula unfenestrated, with 22 mm connector (C) Article number
5	57 mm	1 piece	201402	201432	201462
6	60 mm	1 piece	201403	201433	201463
7	65 mm	1 piece	201404	201434	201464
8	70 mm	1 piece	201405	201435	201465
9	75 mm	1 piece	201406	201436	201466
10	85 mm	1 piece	201407	201437	201467
11	87 mm	1 piece	201408	201438	201468
12	90 mm	1 piece	201409	201439	201469
13	100 mm	1 piece	201410	201440	201470

REPLACEMENT INNER
CANNULAS

Fenestrated replacement inner cannulas
for PRIMEDISTOM® tracheostomy tube
in standard length

(A)



(B)



(C)



SIZE CHART:

Size	Suitable for Primedistom® outer cannula with Length Lb	PU	Inner cannula fenestrated, without connector (A)	Inner cannula fenestrated, with 15 mm connector (B)	Inner cannula fenestrated, with 22 mm connector (C)
			Article number	Article number	Article number
5	57 mm	1 piece	201413	201443	201473
6	60 mm	1 piece	201414	201444	201474
7	65 mm	1 piece	201415	201445	201475
8	70 mm	1 piece	201416	201446	201476
9	75 mm	1 piece	201417	201447	201477
10	85 mm	1 piece	201418	201448	201478
11	87 mm	1 piece	201419	201449	201479
12	90 mm	1 piece	201420	201450	201480
13	100 mm	1 piece	201421	201451	201481



TRACHEOSTOMY TUBES

SILVER tracheostomy tubes are HAND-MADE in GERMANY from high-grade sterling silver (925) of medical quality.

Due to their conical and seamless tube form these tracheostomy tubes enable easy insertion and avoid irritations to the tracheal wall (special versions, e.g. tracheostomy tubes with cylindrical form, are explicitly marked in the catalog).

The material of these tubes has an anti-bacterial impact and their smooth surface is easy to clean. The flange of these cannulas features two lateral eyelets to fix the inserted tracheostomy tube with a neck strap.

The extremely thin walls of these tracheostomy tubes provide a maximum lumen and therefore a large capacity for ventilation, even in small-sized cannulas. The material properties allow the increment of the inner diameter without changing the outer diameter. More air can enter and pass through the cannula. This increased inner volume makes silver cannulas superior to thicker medical-grade polymer cannulas in their intake capacity.


In addition, **SILVER** tracheostomy tubes are characterized by their high durability which allows a longer duration of use. Therefore they are often considered as economical long-term tubes.

This type of tracheostomy tube is always presented in sets with the inner cannula and outer cannula fitting together perfectly. This guarantees safe removal of the inner cannula for cleaning purposes or in case of shortness of breath.

The inner cannula can be fixed or released by moving a small silver lever on the outer cannula.

Customized service is available on request.

TRACHEOSTOMY TUBES




MEASUREMENTS FOR SILVER CONICAL:

Size	ID _{OC}	ID _{IC}	OD	Ln Length along centre line (OC)	Lb Length along outer bend (OC)
3	5,4 mm	4,4 mm	7,5 mm	50 mm	55 mm
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm
14	12,7 mm	11,7 mm	15,8 mm	83 mm	90 mm

ID_{OC} = Inside diameter (clear width) at bottom of outer cannula
ID_{IC} = Inside diameter at bottom of inner cannula
OD = Outside diameter of outer cannula behind the flange
IC = Inner cannula
OC = Outer cannula
Ln = Length along centre line of the outer cannula from neck flange to bottom of the tube
Lb = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit



MEASUREMENTS FOR SILVER CYLINDRICAL:

Size	ID _{oc}	ID _{ic}	OD	Ln Length along centre line (OC)	Lb Length along outer bend (OC)
8	7,2 mm	5,7 mm	8,0 mm	53,6 mm	60 mm
9	8,2 mm	6,4 mm	9,0 mm	58,6 mm	65 mm
10	9,2 mm	7,3 mm	10,0 mm	62,1 mm	70 mm
11	10,2 mm	8,2 mm	11,0 mm	71,2 mm	80 mm
12	11,2 mm	9,1 mm	12,0 mm	80,4 mm	90 mm
13	12,2 mm	10,1 mm	13,0 mm	80,4 mm	90 mm

ID_{oc} = Inside diameter (clear width) at bottom of outer cannula
ID_{ic} = Inside diameter at bottom of inner cannula
OD = Outside diameter of outer cannula behind the flange
IC = Inner cannula
OC = Outer cannula
Ln = Length along centre line of the outer cannula from neck flange to bottom of the tube
Lb = Length along outer bend of the outer cannula from neck flange to bottom of the tube

PU = Packaging unit

Special products like halflexible or flexible **SILVER** tracheostomy tubes are specified with their individual measurements on the relevant product pages.
Further variants not listed in the catalog are available upon request.

TRACHEOSTOMY TUBES

Product features

- Outer cannula
- With inner cannula system
- 2 inner cannulas without connector
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy/ Laryngectomy



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200854
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200855
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200856
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200857
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200858
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200859
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200860
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200861
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200862
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	200863
14	12,7 mm	11,7 mm	15,8 mm	83 mm	90 mm	1 set	200864

INCLUDING:



NECK STRAP



OPTIONAL
with Obturator
(see p. 109)

TRACHEOSTOMY TUBES

Product features

- Outer cannula
- With inner cannula system
- 1 inner cannual without connector
- 1 inner cannula with 15 mm connector
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy/ Laryngectomy



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200869
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200870
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200871
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200872
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200873
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200874
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200875
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200876
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200877

INCLUDING:



NECK STRAP



OPTIONAL ACCESSORIES

matching Obturator
optionally available (see p. 109)

TRACHEOSTOMY TUBES

Product features

- Outer cannula
- With inner cannula system
- 1 inner cannual without connector
- 1 inner cannula with 22 mm connector
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy/ Laryngectomy



SIZE CHART:

Size	ID _{oc}	ID _{ic}	OD	Length		PU	Article number
				Ln	Lb		
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200934
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200935
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200936
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200937
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200938
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200939
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200940
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200941
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200942
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	200943

INCLUDING:



NECK STRAP



OPTIONAL ACCESSORIES

matching Obturator
optionally available (see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- 1 fenestrated inner cannula with detachable silver speaking valve enables phonation (secured against falling off by a small silver chain)
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length		PU	Article number
				Ln	Lb		
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200914
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200915
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200916
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200917
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200918
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200919
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200920
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200921
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200922
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	200923
14	12,7 mm	11,7 mm	15,8 mm	83 mm	90 mm	1 set	200924

INCLUDING:



NECK STRAP



OPTIONAL
with Obturator
(see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- With inner cannula system
- 1 fenestrated inner cannula with detachable silver speaking valve enables phonation (secured against falling off by a small silver chain)
- 1 inner cannula with 15 mm connector
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln	Length Lb	PU	Article number
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200884
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200885
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200886
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200887
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200888
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200889
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200890
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200891
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200892
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	200893

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

matching Obturator
optionally available (see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- With inner cannula system
- 1 fenestrated inner cannula with O₂-port and detachable silver speaking valve which enables phonation (secured against falling off by a small silver chain)
- 1 inner cannula with 15 mm connector
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation
- Hypoxia



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Ln	Length Lb	PU	Article number
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200899
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200900
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200901
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200902
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200903
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200904
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200905
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200906
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200907
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	200908
14	12,7 mm	11,7 mm	15,8 mm	83 mm	90 mm	1 set	200909

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

matching Obturator
optionally available (see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- With inner cannula system
- 1 inner cannula without connector
- 1 fenestrated inner cannula with 22 mm connector enables phonation in conjunction with a speaking valve
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln	Length Lb	PU	Article number
4	6,0 mm	5,0 mm	8,2 mm	50 mm	55 mm	1 set	200834
5	6,7 mm	5,7 mm	9,1 mm	54 mm	60 mm	1 set	200835
6	7,4 mm	6,4 mm	9,8 mm	55 mm	60 mm	1 set	200836
7	8,0 mm	7,0 mm	10,9 mm	60 mm	65 mm	1 set	200837
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	200838
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	200839
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	200840
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	200841
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	200842
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	200843

INCLUDING:



NECK STRAP



OPTIONAL ACCESSORIES

matching Obturator
optionally available (see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- 1 fenestrated inner cannula with detachable round speaking valve enables phonation
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln	Length Lb	PU	Article number
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	202838
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	202839
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	202840
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	202841
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	202842
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	202843

INCLUDING:



NECK STRAP



OPTIONAL
with Obturator
(see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- With inner cannula system
- 1 inner cannula with 15 mm connector
- 1 fenestrated inner cannula with detachable round speaking valve which enables phonation and features an O₂-port
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation
- Hypoxia



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln	Length Lb	PU	Article number
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	202858
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	202859
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	202860
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	202861
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	202862
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	202863

INCLUDING:



NECK STRAP



OPTIONAL
ACCESSORIES

matching Obturator
optionally available (see p. 109)

TRACHEOSTOMY TUBES

Product features

- Sieved outer cannula
- 1 fenestrated inner cannula with detachable square speaking valve which enables phonation
- Conical form facilitates insertion
- Large inner lumen
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy
- Phonation



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Ln	Length Lb	PU	Article number
8	8,7 mm	7,7 mm	11,6 mm	66 mm	70 mm	1 set	202890
9	9,4 mm	8,4 mm	12,3 mm	65 mm	75 mm	1 set	202891
10	10,0 mm	9,0 mm	13,2 mm	74 mm	80 mm	1 set	202892
11	10,7 mm	9,7 mm	13,8 mm	84 mm	90 mm	1 set	202893
12	11,4 mm	10,4 mm	14,6 mm	83 mm	90 mm	1 set	202894
13	12,0 mm	11,0 mm	15,0 mm	83 mm	90 mm	1 set	202895

INCLUDING:



NECK STRAP



OPTIONAL
with Obturator
(see p. 109)

TRACHEOSTOMY TUBES

Product features

- Cylindrical outer cannula
- 1 inner cannula without connector
- 1 Obturator with a smooth atraumatic tip facilitates insertion
- Large inner lumen
- Small neck flange (discrete design)
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), stainless steel, antibacterial

Indication

- Tracheostomy/ Laryngectomy



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Ln	Length Lb	PU	Article number
8	7,2 mm	5,7 mm	8,0 mm	53,6 mm	60 mm	1 set	202902
9	8,2 mm	6,4 mm	9,0 mm	58,6 mm	65 mm	1 set	202903
10	9,2 mm	7,3 mm	10,0 mm	62,1 mm	70 mm	1 set	202904
11	10,2 mm	8,2 mm	11,0 mm	71,2 mm	80 mm	1 set	202905
12	11,2 mm	9,1 mm	12,0 mm	80,4 mm	90 mm	1 set	202906
13	12,2 mm	10,1 mm	13,0 mm	80,4 mm	90 mm	1 set	202907

INCLUDING:



NECK STRAP

TRACHEOSTOMY TUBES

Product features

- Halfflexible outer cannula
- 1 halfflexible inner cannula without connector
- Flexible form adapts perfectly to the individual anatomy of the patient
- Available also with sieving/ fenestration (customization)
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy/ Laryngectomy
- Tracheal stenoses/ Granuloma



SIZE CHART:

Size	ID _{OC}	ID _{IC}	OD	Length Lb	PU	Article number
8	7,4 mm	4,6 mm	11,6 mm	135 mm	1 set	202938
9	8,1 mm	5,3 mm	12,3 mm	140 mm	1 set	202939
10	8,7 mm	5,9 mm	13,2 mm	145 mm	1 set	202940
11	9,4 mm	6,6 mm	13,8 mm	155 mm	1 set	202941
12	10,1 mm	7,6 mm	14,6 mm	155 mm	1 set	202922
13	10,7 mm	7,9 mm	15,0 mm	155 mm	1 set	202943

INCLUDING:



NECK STRAP

TRACHEOSTOMY TUBES

Product features

- Flexible tracheostomy tube
- Flexible form adapts perfectly to the individual anatomy of the patient
- Available also with flexible inner cannula (customization)
- Long duration of use due to high durability of silver material
- Material: sterling silver (925), antibacterial

Indication

- Tracheostomy/ Laryngectomy
- Tracheal stenoses/ Granuloma



SIZE CHART:

Size	ID	OD	Length Lb	PU	Article number
8	4,0 mm	8,0 mm	125 mm	1 set	202958
9	5,0 mm	9,0 mm	130 mm	1 set	202959
10	6,0 mm	10,0 mm	135 mm	1 set	202960
11	7,0 mm	11,0 mm	145 mm	1 set	202961
12	8,0 mm	12,0 mm	155 mm	1 set	202962
13	9,0 mm	13,0 mm	155 mm	1 set	202963

INCLUDING:



NECK STRAP

Product features

- Obturator for conical SILVER tracheostomy tubes
- Facilitates insertion
- With smooth atraumatic tip
- Optionally available in addition to a tracheostomy tube set (Obturator is adapted perfectly to each cannula set)
- Material: stainless steel

Indication

- Tracheostomy/ Laryngectomy



SIZE CHART:

Size (suitable for conical SILVER tracheostomy tube of the same size)	PU	Article number
5	1 set	202805
6	1 set	202806
7	1 set	202807
8	1 set	202808
9	1 set	202809
10	1 set	202810
11	1 set	202811
12	1 set	202812
13	1 set	202813
14	1 set	202814

TRACHEOSTOMY TUBES
FOR NEONATES
AND CHILDREN



TRACHEOSTOMY TUBES

FOR NEONATES

AND CHILDREN

Our tracheostomy tubes for neonates and children are specially tailored to the requirements of little patients. They are made of particulary soft, flexible DEHP-free material and formed atraumaticly. Due to their thin-walled design they ensure high air supply, but at the same time retain their high stability of form. The swivelling 15 mm connector of these tracheostomy tubes enables easy connection of ventilation circuits and minimizes the transmission of tensile and shear forces to the tube.



Swivelling 15 mm connector



The cuffed version have an ultra-thin Polyurethan (PU) cuff



Special recesses on the flange enable maximum flexibility




Different flange designs for neo-natal and pediatric patients

TRACHEOSTOMY TUBES

FOR NEONATES

AND CHILDREN



MEASUREMENTS FOR NEO, PED, PEDL, PEDC, PEDCL:

Size	ID	OD	Ln Length along centre line		
			NEO	PED, PEDC	PEDL, PEDCL
3.0	3,0 mm	4,5 mm	30 mm	39 mm	
3.5	3,5 mm	5,2 mm	32 mm	40 mm	
4.0	4,0 mm	5,9 mm	34 mm	41 mm	
4.5	4,5 mm	6,5 mm	36 mm	42 mm	
5.0	5,0 mm	7,1 mm		44 mm	50 mm
5.5	5,5 mm	7,7 mm		46 mm	52 mm
6.0	6,0 mm	8,3 mm			54 mm
6.5	6,5 mm	9,0 mm			56 mm

ID = Inside diameter at bottom of tube
OD = Outside diameter at bottom of tube
Ln = Length along centre line from neck flange to bottom of the tube

PU = Packaging unit

NEONATAL
TRACHEOSTOMY TUBES
WITHOUT CUFF

Product features

- Neonatal tracheostomy tube made of radio-paque material for position control, with swivelling 15 mm connector
- 1 Obturator for easy insertion
- Without inner cannula system in order to provide a large inner lumen
- Cylindrical form ensures greatest possible inner diameter
- For removal of the swivelling connector and disconnecting of accessories the included disconnecting wedge can be used
- Material: PVC, latex-free, DEHP-free, bisphenol A-free, thermo-sensitive

Indication

- Tracheostomy in neonates
- Radiation/ MRI depending on tumour localisation



NEO - NEONATAL
TRACHEOSTOMY TUBES
WITHOUT CUFF

SIZE CHART:

Size	ID	OD	Length along centre line	PU	Article number
3.0	3,0 mm	4,5 mm	30 mm	1 set	200041
3.5	4,5 mm	5,2 mm	32 mm	1 set	200042
4.0	4,0 mm	5,9 mm	34 mm	1 set	200043
4.5	4,5 mm	6,5 mm	36 mm	1 set	200044

INCLUDING:



NECK STRAP

DISCONNECTING WEDGE



OPTIONAL
ACCESSORIES

PEDIATRIC
TRACHEOSTOMY TUBES
WITHOUT CUFF

Product features

- Pediatric tracheostomy tube made of radiopaque material for position control, with swivelling 15 mm connector
- 1 Obturator for easy insertion
- Without inner cannula system in order to provide a large inner lumen
- Cylindrical form ensures greatest possible inner diameter
- For removal of the swivelling connector and disconnecting of accessories the included disconnecting wedge can be used
- Material: PVC, latex-free, DEHP-free, bisphenol A-free, thermo-sensitive

Indication

- Tracheostomy in children
- Radiation/ MRI depending on tumour localisation



PED - PEDIATRIC
TRACHEOSTOMY TUBE
WITHOUT CUFF

SIZE CHART:

Size	ID	OD	Length along centre line	PU	Article number
3.0	3,0 mm	4,5 mm	39 mm	1 set	200051
3.5	3,5 mm	5,2 mm	40 mm	1 set	200052
4.0	4,0 mm	5,9 mm	41 mm	1 set	200053
4.5	4,5 mm	6,5 mm	42 mm	1 set	200054
5.0	5,0 mm	7,1 mm	44 mm	1 set	200055
5.5	5,5 mm	7,7 mm	46 mm	1 set	200056

INCLUDING:



NECK STRAP

DISCONNECTING WEDGE



OPTIONAL
ACCESSORIES

PEDIATRIC
TRACHEOSTOMY TUBES
WITHOUT CUFF
LONG

Product features

- Long pediatric tracheostomy tube made of radio-paque material for position control, with swivelling 15 mm connector
- 1 Obturator for easy insertion
- Without inner cannula system in order to provide a large inner lumen
- Cylindrical form ensures greatest possible inner diameter
- For removal of the swivelling connector and disconnecting of accessories the included disconnecting wedge can be used
- Material: PVC, latex-free, DEHP-free, bisphenol A-free, thermo-sensitive

Indication

- Tracheostomy in children
- Radiation/ MRI depending on tumour localisation
- Tracheal stenoses/ Granuloma



PEDL - PEDIATRIC
TRACHEOSTOMY TUBE
WITHOUT CUFF, LONG

SIZE CHART:

Size	ID	OD	Length along centre line	PU	Article number
5.0	5,0 mm	7,1 mm	50 mm	1 set	200065
5.5	5,5 mm	7,7 mm	52 mm	1 set	200066
6.0	6,0 mm	8,3 mm	54 mm	1 set	200067
6.5	6,5 mm	9,0 mm	56 mm	1 set	200068

INCLUDING:



NECK STRAP



DISCONNECTING WEDGE



OPTIONAL
ACCESSORIES

PEDIATRIC
TRACHEOSTOMY TUBES
WITH CUFF

Product features

- Pediatric tracheostomy tube made of radiopaque material for position control, with particularly thin-walled low-pressure Polyurethan (PU) cuff and swivelling 15 mm connector
- 1 Obturator for easy insertion
- Without inner cannula system in order to provide a large inner lumen
- Cylindrical form ensures greatest possible inner diameter
- For removal of the swivelling connector and disconnecting of accessories the included disconnecting wedge can be used
- Material: PVC, PU, latex-free, DEHP-free, bisphenol A-free, thermo-sensitive

Indication

- Tracheostomy in children
- Ventilation
- Aspiration



PEDC - PEDIATRIC
TRACHEOSTOMY TUBE WITH CUFF

SIZE CHART:

Size	ID	OD	Length along centre line	PU	Article number
4.0	4,0 mm	5,9 mm	41 mm	1 set	200073
4.5	4,5 mm	6,5 mm	42 mm	1 set	200074
5.0	5,0 mm	7,1 mm	44 mm	1 set	200075
5.5	5,5 mm	7,7 mm	46 mm	1 set	200076

INCLUDING:



OPTIONAL
ACCESSORIES

PEDIATRIC
TRACHEOSTOMY TUBES
WITH CUFF
LONG

Product features

- Long pediatric tracheostomy tube made of radio-paque material for position control, with particularly thin-walled low-pressure Polyurethan (PU) cuff and swivelling 15 mm connector
- 1 Obturator for easy insertion
- Without inner cannula system in order to provide a large inner lumen
- Cylindrical form ensures greatest possible inner diameter
- For removal of the swivelling connector and disconnecting of accessories the included disconnecting wedge can be used
- Material: PVC, PU, latex-free, DEHP-free, bisphenol A-free, thermo-sensitive

Indication

- Tracheostomy in children
- Ventilation
- Aspiration
- Tracheal stenoses/ Granuloma




PEDCL - PEDIATRIC
TRACHEOSTOMY TUBE WITH CUFF
LONG

SIZE CHART:

Size	ID	OD	Length along centre line	PU	Article number
5.0	5,0 mm	7,1 mm	50 mm	1 set	200085
5.5	5,5 mm	7,7 mm	52 mm	1 set	200086
6.0	6,0 mm	8,3 mm	54 mm	1 set	200087
6.5	6,5 mm	9,0 mm	56 mm	1 set	200088

INCLUDING:



OPTIONAL
ACCESSORIES

Product features

- Flexible extension adapter for pediatric tracheostomy tubes
- For ventilated patients it facilitates the connection of a tracheostomy tube to a ventilator
- Wide end (female) can be connected directly to to the 15 mm connector of pediatric tracheostomy tubes
- Narrow end with 15 mm connector (male) on the opposite side can be connected to a ventilator
- For disconnecting the flexible extension adapter from the tracheostomy tube use the disconnecting wedge included with every PRIMED pediatric tracheostomy tube

Indication

- Tracheostomy in children
- Ventilation



FLEXIBLE
EXTENSION ADAPTER
FOR A 15 MM CONNECTOR

	PU	Article number
Flexible Extension Adapter for pediatric tracheostomy tubes (with an adapter for a 15 mm connector), sterile	1 piece	200090

Further accessories for neonates and children you find on page 146 and 157-159.




LARYNGECTOMY TUBES

PRIMEDISILK® laryngectomy tubes are made of pure, medical-grade silicone. Due to high flexibility and the outstanding material properties of silicone, the tubes adjust well to the shape of the trachea and assure high wearing comfort.

An atraumatic rounded tip prevents irritation of the mucous membranes in the trachea and facilitates insertion.

PRIMEDISILK® tubes are recommended for laryngectomized patients whose tracheostoma still has to be stabilized, and for tracheostomized patients, who need a placeholder during the weaning process.

LARYNGECTOMY TUBES



MEASUREMENTS FOR PRIMEDISILK®:

Size	ID	OD	Ln Length along centre line		Lb Length along outer bend	
			short version	standard version	short version	standard version
8	9,5 mm	12,0 mm	37,0 mm	54,0 mm	44,0 mm	62,0 mm
9	10,5 mm	13,5 mm	37,0 mm	54,0 mm	44,0 mm	62,0 mm
10	12,0 mm	15,0 mm	37,0 mm	54,0 mm	44,0 mm	62,0 mm
12	13,5 mm	17,0 mm	37,0 mm	54,0 mm	44,0 mm	62,0 mm

ID = Inside diameter (clear width) at bottom of cannula
OD = Outside diameter of cannula behind the flange
Ln = Length along centre line from neck flange to bottom of the tube
Lb = Length along outer bend from neck flange to bottom of the tube

PU = Packaging unit

LARYNGECTOMY TUBES

Product features

- Low-profile laryngectomy tube without connector
- Prevents shrinkage of tissue
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Atraumatic rounded tip facilitates insertion
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomy/ Tracheostomy
- Radiation/ MRI



PRIMEDISILK®



SHORT VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number
			Ln	Lb		
8	9,5 mm	12,0 mm	54 mm	62 mm	1 set	200614
9	10,5 mm	13,5 mm	54 mm	62 mm	1 set	200615
10	12,0 mm	15,0 mm	54 mm	62 mm	1 set	200616
12	13,5 mm	17,0 mm	54 mm	62 mm	1 set	200617

SHORT VERSION				
Length	Ln	Lb	PU	Article number
37 mm	44 mm	1 set		200660
37 mm	44 mm	1 set		200661
37 mm	44 mm	1 set		200662
37 mm	44 mm	1 set		200663

INCLUDING:



NECK STRAP

LARYNGECTOMY TUBES

Product features

- Laryngectomy tube with 22 mm connector
- Prevents shrinkage of tissue
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Atraumatic rounded tip facilitates insertion
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomy/ Tracheostomy
- Radiation/ MRI



PRIMEDISILK-VENT®




SHORT VERSION


SIZE CHART:

Size	ID	OD	Length		PU	Article number	SHORT VERSION			
			Ln	Lb			Length Ln	Lb	PU	Article number
8	9,5 mm	12,0 mm	54 mm	62 mm	1 set	200628	37 mm	44 mm	1 set	200668
9	10,5 mm	13,5 mm	54 mm	62 mm	1 set	200629	37 mm	44 mm	1 set	200669
10	12,0 mm	15,0 mm	54 mm	62 mm	1 set	200630	37 mm	44 mm	1 set	200670
12	13,5 mm	17,0 mm	54 mm	62 mm	1 set	200631	37 mm	44 mm	1 set	200671

INCLUDING:



NECK STRAP



5 X PRIMEDIAIR-PHON
ULTRAFLOW



OPTIONAL
ACCESSORIES

LARYNGECTOMY TUBES

Product features

- Sieved laryngectomy tube with 22 mm connector
- Sieving enables phonation (for laryngectomized patients with normal position of voice prosthesis)
- Prevents shrinkage of tissue
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Atraumatic rounded tip facilitates insertion
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomized patients (with voice prosthesis)/
Tracheostomized patients with intact larynx
- Phonation
- Radiation/ MRI



PRIMEDISILK-VENT®



SHORT VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number	SHORT VERSION			
			Ln	Lb			Length Ln	Lb	PU	Article number
8	9,5 mm	12,0 mm	54 mm	62 mm	1 set	200624	37 mm	44 mm	1 set	200664
9	10,5 mm	13,5 mm	54 mm	62 mm	1 set	200625	37 mm	44 mm	1 set	200665
10	12,0 mm	15,0 mm	54 mm	62 mm	1 set	200626	37 mm	44 mm	1 set	200666
12	13,5 mm	17,0 mm	54 mm	62 mm	1 set	200627	37 mm	44 mm	1 set	200667

INCLUDING:

NECK STRAP

5 X PRIMEDIAIR-PHON
ULTRAFLW



OPTIONAL
ACCESSORIES

LARYNGECTOMY TUBES

Product features

- Laryngectomy tube with 22 mm connector without flange
- To be used in conjunction with adhesive base-plates (PRIMEDISTRIP® stoma patches)
- Enables application without neck strap
- Prevents shrinkage of tissue
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Atraumatic rounded tip facilitates insertion
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomy/ Tracheostomy
- Radiation/ MRI



PRIMEDISILK-CLIP®



SHORT VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number	SHORT VERSION			
			Ln	Lb			Ln	Lb	PU	Article number
8	9,5 mm	12,0 mm	54 mm	62 mm	1 set	200684	37 mm	44 mm	1 set	200692
9	10,5 mm	13,5 mm	54 mm	62 mm	1 set	200685	37 mm	44 mm	1 set	200693
10	12,0 mm	15,0 mm	54 mm	62 mm	1 set	200686	37 mm	44 mm	1 set	200694
12	13,5 mm	17,0 mm	54 mm	62 mm	1 set	200687	37 mm	44 mm	1 set	200695



OPTIONAL ACCESSORIES

LARYNGECTOMY TUBES

Product features

- Sieved laryngectomy tube with 22 mm connector without flange
- Sieving enables phonation (for laryngectomized patients with normal position of voice prosthesis)
- To be used in conjunction with adhesive base-plates (PRIMEDISTRIP® stoma patches)
- Enables application without neck strap and an improved sealing of the tracheostoma during phonation
- Prevents shrinkage of tissue
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Atraumatic rounded tip facilitates insertion
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomized patients (with voice prosthesis)/ Tracheostomized patients with intact larynx
- Phonation
- Radiation/ MRI



PRIMEDISILK-CLIP®
WITH SIEVING



SHORT VERSION

SIZE CHART:

Size	ID	OD	Length		PU	Article number
			Ln	Lb		
8	9,5 mm	12,0 mm	54 mm	62 mm	1 set	200680
9	10,5 mm	13,5 mm	54 mm	62 mm	1 set	200681
10	12,0 mm	15,0 mm	54 mm	62 mm	1 set	200682
12	13,5 mm	17,0 mm	54 mm	62 mm	1 set	200683

SHORT VERSION			
Length	Ln	Lb	PU
	37 mm	44 mm	1 set
	37 mm	44 mm	1 set
	37 mm	44 mm	1 set
	37 mm	44 mm	1 set



OPTIONAL
ACCESSORIES

TRACHEOSTOMY PLACEHOLDERS

PRIMEDISK® tracheostomy placeholders provide a secondary airway by gently holding the tracheostoma open. They provide access during a bronchial toilet and are used for decannulation.

Made of soft, medical-grade silicone they ensure a very high wearing comfort and good tissue compatibility.

We offer different types of **PRIMEDISK®** tracheostomy placeholders which should be selected according to the individual anatomical features of each patient:

TRACHEO-SAFE T-TUBUS STOMA-BUTTON



PU = Packaging unit
ID = Inside diameter
OD = Outside diameter

Product features

- Tracheostomy placeholder following a decannulation and prior to a tracheostomy seal
- The extended semitubular form of the tracheal shank ideally adapts to the shape of the trachea
- Variably adjustable skin-colored flange ensures secure positioning in the trachea
- Optional closure plug enables respiration via the upper respiratory tract (can be opened with ease at any time for suctioning of secretion)
- Prevents shrinkage of tissue
- Enables patient to cough up the stimulus secretion very fast
- Provides an artificial airway and access for pulmonary hygiene
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Material: silicone, latex-free

Indication

- Tracheostomy
- Patients with stable breathing who primarily breath through the upper airways
- Weaning
- Radiation/ MRI

SIZE CHART:

Length	ID	OD	PU	Article number
45 mm	8 mm	14,0 mm	1 piece	20895
70 mm	8 mm	13,8 mm	1 piece	20898



TRACHEO-SAFE


T-TUBUS

Product features

- Tracheostomy placeholder for the stabilization of the trachea
- Maintains an adequate airway from larynx to mouth and nose
- Height-adjustable distancing ring ensures optimum fit and protects against dislocation
- Can be closed off, if necessary, with the sealing plug (removing the plug allows suction of tracheal secretions)
- Atraumatically shaped ends ensure easy insertion
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Material: silicone, latex-free

Indication

- Tracheostomy
- Tracheal and subglottic stenoses and related reconstructive measures
- Tracheomalacia/ unstable trachea
- Radiation/ MRI



T-TUBUS

SIZE CHART:

Size	A	B	C	D (OD)	E (OD)	PU	Article number
8	63 mm	49 mm	45 mm	6,5 mm	8 mm	1 piece	200640
9	75 mm	50 mm	63 mm	8 mm	9 mm	1 piece	200641
10	75 mm	50 mm	63 mm	8 mm	10 mm	1 piece	200642
11	75 mm	50 mm	63 mm	9 mm	11 mm	1 piece	200643
12	75 mm	75 mm	48 mm	11 mm	12 mm	1 piece	200644
13	100 mm	50 mm	63 mm	11 mm	13 mm	1 piece	200645
14	100 mm	50 mm	63 mm	11 mm	14 mm	1 piece	200646

Product features

- Tracheostomy placeholder for the stabilization of the tracheostoma
- Flat unobtrusive design with ring-shaped shield
- Integrated bulge at the proximal end ensures that it is held securely in the trachea
- Prevents shrinkage of tissue
- Short, atraumatic form can be introduced with ease and reduces the associated risk of injury
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Can be closed with a finger to enable phonation
- Material: silicone, latex-free

Indication

- Laryngectomized patients (with voice prosthesis)/ Tracheostomized patients with intact larynx
- Weaning
- Radiation/ MRI



STOMA BUTTON
TYPE I

SIZE CHART:

Size	ID	OD	Length	PU	Article number
6	8,0 mm	9,5 mm	12,5 mm	1 piece	200602
8	9,5 mm	12,0 mm	13,5 mm	1 piece	200603
9	10,5 mm	13,5 mm	14,5 mm	1 piece	200604
10	12,0 mm	15,0 mm	15,5 mm	1 piece	200605
10 (short)	12,0 mm	15,0 mm	6,5 mm	1 piece	200606
12	13,5 mm	17,0 mm	14,5 mm	1 piece	200607

Product features

- Tracheostomy placeholder for the stabilization of the tracheostoma
- The shield is a funnel-shaped ring which ensures effective sealing against the tracheostoma
- Integrated bulge at the proximal end ensures that it is held securely in the trachea
- Available with and without eyelets (for additional fixation with a neck strap)
- Prevents shrinkage of tissue
- Short, atraumatic form can be introduced with ease and reduces the associated risk of injury
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- Can be closed with a finger to enable phonation
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomized patients (with voice prosthesis)/ Tracheostomized patients with intact larynx
- Weaning
- Radiation/ MRI



STOMA BUTTON
TYPE II



SIZE CHART:

Size	ID	OD	Length	PU	WITHOUT EYELETS (A)	WITH EYELETS (B)
					Article number	Article number
8	8,0 mm	11,0 mm	15,0 mm	1 piece	200555	200545
			22,0 mm	1 piece	200550	200540
9	9,0 mm	12,0 mm	15,0 mm	1 piece	200556	200546
			22,0 mm	1 piece	200551	200541
10	10,0 mm	13,0 mm	15,0 mm	1 piece	200557	200547
			22,0 mm	1 piece	200552	200542
11	11,0 mm	14,0 mm	15,0 mm	1 piece	200558	200548
			22,0 mm	1 piece	200553	200543
12	12,0 mm	15,0 mm	15,0 mm	1 piece	200559	200549
			22,0 mm	1 piece	200554	200544

Product features

- Tracheostomy placeholder for the stabilization of the tracheostoma
- The shield is a funnel-shaped housing which ensures effective sealing against the tracheostoma and features an integrated 22 mm connector
- Integrated bulge at the proximal end ensures that it is held securely in the trachea
- Available with and without eyelets (for additional fixation with a neck strap)
- Prevents shrinkage of tissue
- Short, atraumatic form can be introduced with ease and reduces the associated risk of injury
- High wearing comfort due to flexible and soft material
- Ensures good tissue compatibility
- In order to enable phonation it can either be closed with a finger or a speaking valve can be attached
- Available in two lengths
- Material: silicone, latex-free

Indication

- Laryngectomized patients (with voice prosthesis)/ Tracheostomized patients with intact larynx
- Weaning
- Radiation/ MRI

SIZE CHART:

Size	ID	OD	Length	PU	WITHOUT EYELETS (A)	WITH EYELETS (B)
					Article number	Article number
8	8,0 mm	11,0 mm	15,0 mm	1 piece	200575	200565
			22,0 mm	1 piece	200570	200560
9	9,0 mm	12,0 mm	15,0 mm	1 piece	200576	200566
			22,0 mm	1 piece	200571	200561
10	10,0 mm	13,0 mm	15,0 mm	1 piece	200577	200567
			22,0 mm	1 piece	200572	200562
11	11,0 mm	14,0 mm	15,0 mm	1 piece	200578	200568
			22,0 mm	1 piece	200573	200563
12	12,0 mm	15,0 mm	15,0 mm	1 piece	200579	200569
			22,0 mm	1 piece	200574	200564



STOMA BUTTON
TYPE III



OPTIONAL
ACCESSORIES

PRIMEDISTRIP® is a retaining ring in combination with a self-adhesive patch as a carrier plate.

PRIMEDISTRIP® is suitable for patients with a tracheostoma who, although no longer needing a tracheostomy tube, still require aids such as heat and moisture exchangers (HME) or speaking valves.

It can be secured easily and rapidly, thanks to the skin-friendly adhesive surface. It provides secure protection against the escape of secretion and air. Irritation of the skin is avoided and an optimum seal achieved as a result.

The retaining ring of a **PRIMEDISTRIP®** has an inner diameter of 22 mm and enables attachment of speaking valves (e.g. PRIMEDI RESISTANCE VALVE and PRIMEDI LOW-RESISTANCE VALVE) or heat and moisture exchangers (e.g. PRIMEDIAIR® KAS, PRIMEDIAIR® PHON or PRIMEDIAIR® LATERAL [Type C and D]).

PRIMEDISTRIP® stoma patches can be selected in different variants (Basic, Basic Plus, Comfort, Comfort Plus, Sensitive) to suit the individual requirements of each patient. These are available in different shapes (oval, round), materials and adhesive strengths.




Product features


- Self-adhesive stoma patch with a 22 mm retaining ring for sealing the tracheostoma
- Enables attachment of heat and moisture exchangers or speaking vales
- Transparent
- Easy to remove
- Available in round and oval shape
- Average utilization period: 24 hours

Indication

- Tracheostomy/ Laryngectomy
- Suitable for patients with a normal stoma and normal skin



BASIC



	PU	Shape	Article number
A	15 pieces	round	200360
B	15 pieces	oval	200362




OPTIONAL ACCESSORIES

Product features


- Self-adhesive stoma patch with a 22 mm retaining ring for sealing the tracheostoma
- Enables attachment of heat and moisture exchangers or speaking vales
- Transparent
- Highly flexible
- Very strong adhesive properties
- Available in round and oval shape
- Average utilization period: 48 hours

Indication

- Tracheostomy/ Laryngectomy
- Suitable for patients with an irregular or deep stoma



BASIC PLUS



	PU	Shape	Article number
A	15 pieces	round	200364
B	15 pieces	oval	200365



OPTIONAL ACCESSORIES

Product features


- Self-adhesive stoma patch with a 22 mm retaining ring for sealing the tracheostoma
- Enables attachment of heat and moisture exchangers or speaking vales
- Skin-friendly hydrocolloid material absorbs moisture, thereby ensuring high wearing comfort
- Skin-colored
- Combines high flexibility with strong adhesive properties
- Available in round and oval shape
- Average utilization period: 36 hours

Indication

- Tracheostomy/ Laryngectomy
- Suitable for patients with normal and sensitive skin



COMFORT



	PU	Shape	Article number
A	15 pieces	round	200361
B	15 pieces	oval	200363




OPTIONAL ACCESSORIES

Product features


- Self-adhesive stoma patch with a 22 mm retaining ring for sealing the tracheostoma
- Enables attachment of heat and moisture exchangers or speaking vales
- Soft texture ensures high wearing comfort
- Skin-colored
- Highly flexible
- Strong adhesive properties
- Available in round and oval shape
- Average utilization period: 48 hours

Indication

- Tracheostomy/ Laryngectomy
- Suitable for patients with normal and sensitive skin



COMFORT PLUS



	PU	Shape	Article number
A	15 pieces	round	200366
B	15 pieces	oval	200367




OPTIONAL ACCESSORIES

Product features


- Self-adhesive stoma patch with a 22 mm retaining ring for sealing the tracheostoma
- Made of special hypoallergenic material that avoids skin irritations and ensures high wearing comfort
- White
- Highly flexible
- Available in round and oval shape
- Average utilization period: 48 hours

Indication

- Tracheostomy/ Laryngectomy
- Suitable for patients with normal and highly sensitive skin



SENSITIVE



	PU	Shape	Article number
A	15 pieces	round	200369
B	15 pieces	oval	200370




Product features

- Sealing patch of skin-coloured, biocompatible silicone used for sealing and padding between the tracheostoma and a tracheostomy tube
- Skin-friendly adhesive material enables rapid fixing without adhesives
- Highly flexible
- Natural adaptation to the musculature and tissue allows for a new quality of freedom of movement
- Due to its plano-convex shape the patch even adapts to irregular tracheostomata and seals these effectively
- Improves phonation through the high degree of sealing (making it also suitable for laryngectomized patients with voice prostheses)
- Stabilizes the seat of the tracheostomy tube
- Reduces the emergence of air and secretion, thus improves secretion management
- Facilitates the setting of uniform ventilation parameters in the case of ventilated patients
- Skin-colored
- Available in three sizes, for adults and children

Indication

- Tracheostomy/ Laryngectomy
- Suitable for patients with normal and sensitive skin
- Suitable for patients with an irregular stoma



	PU	Shape	Article number
PRIMED TRACHEO SEAL® for adults, round (A)	3 pieces	round	200236
PRIMED TRACHEO SEAL® for adults, oval (B)	3 pieces	oval	200237
PRIMED TRACHEO SEAL® for children, round (C)	3 pieces	round	200238

Accessories for pulmonary rehabilitation are specially designed for the needs of patients both in hospital and rehabilitation care and receiving nursing care at home.

Patients are either completely or partially incapable of breathing via their larynx, mouth and nose following a laryngectomy or tracheostomy. This is also accompanied by a loss of function of the upper respiratory tract. The accessories specially developed for laryngectomized and tracheostomized patients take on these functions, achieving a considerable improvement in the quality of life.

In addition to **heat and moisture exchangers (HME)** which contribute to pulmonary rehabilitation through the humidification, heating and filtration of respiratory air, we also provide a large selection of **speaking valves** which are suitable for both laryngectomized patients with voice prosthesis and tracheostomized patients with an intact larynx which enable phonation. We also provide **speaking valves with a filter function** and **respiration valves** (for ventilated tracheostomized patients).



PU = Packaging unit
I = Inspiration
Ex = Expiration

Product features

- Effective filtration, humidification and heating of inhaled air (ensures a well-functioning mucous cilia system for infection prevention)
- Low breathing resistance thanks to large surface
- With foam filter
- Standardized 15 mm opening for attachment to a tracheostomy tube or endotracheal tube with 15 mm connection
- Integrated oxygen supply port
- PRIMEDIAIR® VO2 features additionally a self-opening safety/ suction port in the lid which can be used for suctioning and in case of coughing or blockage acts as a relief valve

Indication

- Spontaneously breathing intubated, tracheostomized and laryngectomized patients



PRIMEDIAIR®-VO1



PRIMEDIAIR®-VO2

	PU	Connector	Article number
PRIMEDIAIR®-VO1	30 pieces	15 mm	21066
PRIMEDIAIR®-VO2	30 pieces	15 mm	21063

PRIMEDIAIR®-VO1/ PRIMEDIAIR®-VO2		
Resistance at		
30 l/min	0,08 hPa (I);	0,09 hPa (Ex)
60 l/min	0,13 hPa (I);	0,14 hPa (Ex)
90 l/min	0,16 hPa (I);	0,17 hPa (Ex)
Moisture loss		
Vt = 250 ml	8,5 mg/l	
Vt = 500 ml	10,6 mg/l	
Vt = 750 ml	13,0 mg/l	
Vt = 1000 ml	17,0 mg/l	

Product features

- Effective filtration, humidification and heating of inhaled air (ensures a well-functioning mucous cilia system for infection prevention)
- Low breathing resistance thanks to large surface
- With 15 mm connector for attachment to a tracheostomy tube or endotracheal tube with 15 mm connection
- Integrated oxygen supply port
- Integrated self-opening safety/ suction port in the lid which can be used for suctioning and in case of coughing or blockage acts as a relief valve
- Available with three different filter media:

- Type A with hygroscopic foam filter
- Type B with hygroscopic UltraFlow foam filter with extra low breathing resistance (especially suitably during physical exertion)
- Type C with hygroscopic microwell paper

Indication

- Spontaneously breathing intubated, tracheostomized and laryngectomized patients



PRIMEDIAIR®-DUAL

	PU	Connector	Article number
WITH FOAM FILTER (A)	30 pieces	15 mm	21065
WITH ULTRAFLOW FOAM FILTER (B)	30 pieces	15 mm	200352
WITH PAPER FILTER (C)	30 pieces	15 mm	20497

	WITH FOAM FILTER (A)	WITH ULTRAFLOW-FOAM FILTER (B)	WITH PAPER FILTER (C)
Resistance at			
30 l/min	0,30 hPa (I)	0,06 hPa (I)	0,32 hPa (I)
	0,34 hPa (Ex)	0,08 hPa (Ex)	0,28 hPa (Ex)
60 l/min	0,69 hPa (I)	0,20 hPa (I)	0,82 hPa (I)
	0,78 hPa (Ex)	0,20 hPa (Ex)	0,68 hPa (Ex)
90 l/min	1,12 hPa (I)	0,40 hPa (I)	1,42 hPa (I)
	1,28 hPa (Ex)	0,35 hPa (Ex)	1,14 hPa (Ex)
Moisture loss			
Vt = 250 ml	9,7 mg/l	9,41 mg/l	12,10 mg/l
Vt = 500 ml		11,02 mg/l	14,68 mg/l
Vt = 750 ml		12,82 mg/l	15,91 mg/l
Vt = 1000 ml		13,30 mg/l	17,36 mg/l

Product features

- Effective filtration, humidification and heating of inhaled air (ensures a well-functioning mucous cilia system for infection prevention)
- Special design that allows the exhaled air to pass through the side openings, thereby avoiding exposure of exhaled air to others
- Ensures an optimum air supply, even when covered with a cloth or tracheostoma protection
- With foam filter
- Integrated oxygen supply port
- Available in two different connection versions, with 15 mm connector (A, B) or 22 mm connector (C, D)
- In addition, type B and D are equipped with an integrated self-opening safety/ suction port in the lid which can be used for suctioning and in case of coughing or blockage acts as a relief valve

Indication

- Spontaneously breathing tracheostomized and laryngectomized patients

	PU	Connector	Safety/ suction port	Article number
A	30 pieces	15 mm	--	20492
B	30 pieces	15 mm	✓	20493
C	30 pieces	22 mm	--	20494
D	30 pieces	22 mm	✓	20495

PRIMEDIAIR®-LATERAL (A,B,C,D)		
Resistance at		
30 l/min	0,6 hPa (I);	1,1 hPa (Ex)
60 l/min	2,5 hPa (I);	3,3 hPa (Ex)
90 l/min	5,6 hPa (I);	6,2 hPa (Ex)
Moisture loss		
Vt = 250 ml	11,8 mg/l	
Vt = 500 ml	15,6 mg/l	
Vt = 750 ml	16,8 mg/l	
Vt = 1000 ml	18,1 mg/l	



PRIMEDIAIR®-LATERAL



Product features

- Effective filtration, humidification and heating of inhaled air (ensures a well-functioning mucous cilia system for infection prevention)
- Through several round recesses in the lid air can be inhaled and exhaled directly from the front
- With 15 mm connector for attachment to a tracheostomy tube with 15 mm connection
- With foam filter
- PRIMEDIAIR-O₂ features additionally an integrated oxygen supply port

Indication

- Spontaneously breathing tracheostomized and laryngectomized patients

	PU	Connector	Article number
PRIMEDIAIR®-Kan	30 pieces	15 mm	20485
PRIMEDIAIR®-O2	30 pieces	15 mm	20486

PRIMEDIAIR®-Kan/ PRIMEDIAIR®-O2	
Resistance at	
30 l/min	0,6 cm H ₂ O (hPa)
60 l/min	2,1 cm H ₂ O (hPa)
90 l/min	4,3 cm H ₂ O (hPa)
Moisture loss	
Vt = 250 ml	10,3 mg/l
Vt = 500 ml	13,5 mg/l
Vt = 750 ml	16,3 mg/l
Vt = 1000 ml	18,6 mg/l



PRIMEDIAIR®-KAN



PRIMEDIAIR®-O2

Product features

- Effective filtration, humidification and heating of inhaled air (ensures a well-functioning mucous cilia system for infection prevention)
- Particularly flat and unobtrusive design
- With foam filter
- With 22 mm connector for attachment to a tracheostomy or laryngectomy tube with 22 mm connection, a stoma button (Type III) or a stoma patch (PRIMEDISTRIP®)

Indication

- Spontaneously breathing tracheostomized and laryngectomized patients



PRIMEDIAIR®-KAS

	PU	Connector	Article number
PRIMEDIAIR®-Kas	30 pieces	22 mm	20482

Product features

- Effective filtration, humidification and heating of inhaled air (ensures a well-functioning mucous cilia system for infection prevention)
- Specially designed for the tidal volume of neonates and children
- With 15 mm connector for attachment to a tracheostomy tube with 15 mm connection
- Available in two different filter media:

Typ A with hygroscopic blue foam filter
Typ B with hygroscopic white paper filter

Indication

- Spontaneously breathing tracheostomized neonates and children



PRIMEDIAIR®-NATAL



	PU	Connector	Article number
A With foam filter	30 pieces	15 mm	21062
B With paper filter	30 pieces	15 mm	20496

	PRIMEDIAIR®- Natal with foam filter (A)	PRIMEDIAIR®- Natal with paper filter (B)
Restistance at		
5 l/min	0,14 hPa (I) 0,14-0,16 hPa (Ex)	0,3 hPa (I) 0,3 hPa (Ex)
Moisture loss		
Vt = 30 ml	10,2 mg/l	9,7 mg/l

Product features

- Enables attachment of an HME to an oxygen unit
- Funnel connector can be connected to the oxygen unit; Special connector with tube piece can be attached to the oxygen supply port of the HME
- Internal profile prevents kinking
- Ch.14; Length 200 cm
- Individually packaged, sterile
- Material: PVC, DEHP-free, latex-free

Indication

- Patients who need an external supply of oxygen/
Hypoxia



OXYGEN
CONNECTING TUBE

	PU	Connection	Article number
Oxygen connecting tube for HME, sterile	50 pieces	Funnel - Tube piece	20796

Speaking valves can be divided into two different groups. There are valves for tracheostomized patients with preserved larynx and those for laryngectomized patients with voice prosthesis. There are technical and functional differences between both types of speaking valves.

Speaking valves for tracheostomized patients with preserved larynx:

Tracheostomized patients can use the natural respiratory tract via the mouth and nose, breathing in through the speaking valve and exhaling via the upper respiratory tract. As the speaking valve closes following inspiration, the expiration air escapes over the available larynx, thus enabling phonation. Speaking valves suitable for this purpose include the PRIMEDI-PHON®, PRIMEDI SPEAKING VALVE, PRIMEDI-PHON PLUS®, PRIMEDI-PHON VENT® and the SILICONE SPEAKING VALVE.

These speaking valves for tracheostomized patients must not be used on laryngectomees!

Speaking valves for laryngectomized patients with voice prosthesis:

Laryngectomized patients with a voice prosthesis can cultivate their own compensatory voice with the aid of a special speaking valve. As these patients breathe exclusively via the tracheostoma, the speaking valve needs to remain open during the respiration process to ensure an unobstructed exchange of air. For this reason, speaking valves for laryngectomized patients are either purposely closed with a finger to generate speech (PRIMEDIAIR®-PHON) or function with a special valve membrane which, due to the increased exhalation pressure during phonation, diverts the expiration air over the voice prosthesis, but otherwise remains open (PRIMEDI RESISTANCE VALVES).

These speaking valves can also be used by tracheostomized patients with an intact larynx.



PRIMEDI-PHON®

Product features

- Handsfree speaking valve
- Skin-colored
- With 22 mm connector for attachment to sieved/ fenestrated PRIMEDISILK®, PRIMEDISTOM® or SILVER tubes with 22 mm connection

Indication

- Spontaneously breathing tracheostomized patients
(The use on laryngectomized patients is contraindicated!)

	PU	Connector	Article number
PRIMEDI-PHON®	1 piece	22 mm	200358



PRIMEDI-PHON®



PRIMEDI SPEAKING VALVE

PRIMEDI SPEAKING VALVE

Product features

- Handsfree speaking valve
- White design
- With 15 mm connector for attachment to sieved/ fenestrated tracheostomy tubes with 15 mm connection (suitable for all tube sizes)

Indication

- Spontaneously breathing tracheostomized patients
(The use on laryngectomized patients is contraindicated!)

	PU	Connector	Article number
PRIMEDI SPEAKING VALVE	1 piece	15 mm	20489

Product features

- Handsfree speaking valve with closed basic position of the silicone membrane
- Available in two different colour schemes:

transparent: suitable for the use in the outpatient sector and for home care, due to its unobtrusive colour scheme

blue-transparent: designed for use with inpatients, due to its conspicuous colouring
- Available with or without integrated oxygen supply port
- Suitable for the following:
Sieved/ fenestrated tracheostomy tubes
Tracheostomy tubes with completely unblocked cuff and tracheostomy tubes without cuff

(IMPORTANT: An adequate airflow around the tracheostomy tube must be ensured!)
- The filter cartridge contained in the scope of delivery can be optionally connected to the PRIMEDI-PHON PLUS® speaking valve for filtering inspiration air (also available separately)

How does it work?

- A special feature of this speaking valve is the closed basic position of the silicone membrane
- The membrane only opens briefly during inspiration and reacts to even very weak inhalation
- Afterwards it automatically returns to the tightly sealed basic position without exhalation pressure (expiration pressure) being required (an extremely short shutting time); phonation therefore does not require any build-up of expiration pressure



**PRIMEDI-PHON PLUS®
(TRANSPARENT)**



SPEAKING VALVES
FOR TRACHEOSTOMIZED
PATIENTS

- The impermeability of the PRIMEDI-PHON PLUS® speaking valve creates an air column in the tube during exhalation (in sieved tubes above the sieving) which prevents secretion from penetrating and, consequently, any relocation of the tracheostomy tube
- The expiration air is channelled completely via the larynx and upper respiratory tract
- This enables the formation of speech (phonation) and expectoration of secretions

Indication

- Spontaneously breathing tracheostomized patients
(The use on laryngectomized patients is contraindicated!)



PRIMEDI-PHON PLUS®
(BLUE-TRANSPARENT)



	PU	O ₂ -port	Connector	Article number
PRIMEDI-PHON® PLUS (transparent)	1 piece	no	15 mm/ 22 mm	200354
PRIMEDI-PHON® PLUS (transparent, with O ₂ -port)	1 piece	yes	15 mm/ 22 mm	200355
PRIMEDI-PHON® PLUS (blue-transparent)	1 piece	no	15 mm/ 22 mm	200356
PRIMEDI-PHON® PLUS (blue-transparent, with O ₂ -port)	1 piece	yes	15 mm/ 22 mm	200357
FILTER CARTRIDGE 15/22 for PRIMEDI-PHON® PLUS	30 pieces	no	15 mm	200351

INCLUDING:



SPEAKING VALVES
FOR TRACHEOSTOMIZED
PATIENTS

Product features

- Handsfree speaking valve
- Made of medical grade silicone
- Operates noiselessly
- Exclusivly compatible with the fenestrated PRIMEDISTOM® inner cannula without connector size 5 to 13

Indication

- Spontaneously breathing tracheostomized patients
(The use on laryngectomized patients is contraindicated!)



SILICONE
SPEAKING VALVE

	PU	Connector	Article number
Silicone speaking valve	1 piece	compatible with PRIMEDISTOM® inner cannula without connector, size 5-13	200350

Product features

- Handsfree speaking valve with closed basic position of the silicone membrane and connection for attachment to a ventilator
- Suitable for the following:
Sieved/ fenestrated tracheostomy tubes
Tracheostomy tubes with completely unblocked cuff and tracheostomy tubes without cuff
(IMPORTANT: An adequate airflow around the tracheostomy tube must be ensured!)
- Can be attached to a ventilator via swivelling 15 mm connector (For setting of ventilation parameters observe instructions for use!)
- Inhalation takes place via ventilator; at the end of inspiration the silicone membrane automatically returns to the tightly sealed basic position and the expiration air is channelled completely via the larynx and upper respiratory tract, thus enabling phonation

Indication

- Tracheostomized, ventilated patients with intact larynx (adults, children and neonates)
- Spontaneously breathing tracheostomized patients (**The use on laryngectomized patients is contraindicated!**)
- Neuromuscular disorders
- Ventilation
- Weaning



PRIMEDI-PHON VENT® (PPV)



	PU	Connector	Article number
PRIMEDI-PHON® VENT	1 piece	15 mm/ 22 mm (near the patient) - swivelling 15 mm (near the ventilator)	200320

INCLUDING:



15 MM TO 22 MM
ADAPTER

**SPEAKING VALVES
FOR LARYNGECTOMIZED &
TRACHEOSTOMIZED
PATIENTS**

Product features

- Handsfree speaking valve with silicone membrane
- Infinitely adjustable (via the integrated rotation mechanism) to the individual exhalation pressure of the patient
- Available in two different versions:
 - A PRIMEDI-RESISTANCE-VALVE with normal resistance (blue silicone membrane)
 - B PRIMEDI-LOW-RESISTANCE-VALVE with low resistance (white silicone membrane); here only relatively little pressure is required to close it
- Air can pass unobstructed through the speaking valve during respiration. Only the increased exhalation pressure during expiration leads to closing of the silicone membrane, and exhalation air is conducted through the voice prosthesis or the larynx to generate speech
- With 22 mm connector for attachment to a tracheostomy or laryngectomy tube with 22 mm connection, a stoma button (Type III) or a stoma patch (PRIMEDISTRIP®)

Indication

- Laryngectomized patients (with voice prosthesis)/
Tracheostomized patients with intact larynx



PRIMEDI-RESISTANCE-VALVES



	PU	Connector	Article number
PRIMEDI-RESISTANCE-VALVE (A)	1 piece	22 mm	20481
PRIMEDI-LOW-RESISTANCE-VALVE (B)	1 piece	22 mm	20480

SPEAKING VALVES WITH FILTER FUNCTION
FOR LARYNGECTOMIZED &
TRACHEOSTOMIZED
PATIENTS

Product features

- Speaking valve with filter function
- Low profile
- Discrete design: skin-coloured or transparent
- Filtration, humidification and heating of inhaled air thanks to the integrated special foam (resulting in a reduction of secretion, coughing and incrustations)
- Enables phonation merely by closing it with a finger
- With 22 mm connector for attachment to a tracheostomy or laryngectomy tube with 22 mm connection, a stoma button (Type III) or a stoma patch (PRIMEDISTRIP®)
- Available in two different filter media:
PRIMEDIAIR®-PHON with normal foam filter
PRIMEDIAIR®-PHON ULTRAFLOW with Ultra-Flow foam filter with a particularly low resistance to inspiration
- PRIMEDIAIR®-PHON ULTRAFLOW facilitates respiration and is particularly recommended for strenuous physical activity

Indication

- Laryngectomized patients (with voice prosthesis)/
Tracheostomized patients with intact larynx



PRIMEDIAIR®-PHON



PRIMEDIAIR®-PHON
ULTRAFLOW

	PU	Connector	Article number
PRIMEDIAIR® -PHON	30 pieces	22 mm	20483
PRIMEDIAIR® -PHON ULTRAFLOW	30 pieces	22 mm	200359

Product features

- BACTERIA and VIRUS FILTERS for protecting the patient as well as the ventilation unit during ventilation
- Verified retention effect of 99,999% for bacteria and viruses
- Minimize the risk of cross-infections
- Feature a built-in CO₂- port
- Designed for adults with a tidal volume of 300-1500 ml
- Available with and without integrated heat and moisture exchanger (HME)
- Filter with HME additionally humidifies and warms the respiratory gases
- To be positioned between an endotracheal tube or tracheostomy tube and the Y-piece of a ventilation system

Indication

- Ventilation

(Do not use in combination with volatile oils and inhalation devices!)



WITHOUT HME



WITH HME

	WITHOUT HME	WITH HME
Article number	21218	21219
Connector	22 mm F/ 15 mm M	22 mm F/ 15 mm M
	22 mm M/15 mm F	22 mm M/15 mm F
PU	50 pieces	50 pieces
Compressible volume	45 ml	45 ml
Resistance at		
30 l/min	1,1 hPa (I)	1,4 hPa (I)
	1,2 hPa (Ex)	1,3 hPa (Ex)
60 l/min	2,6 hPa (I)	3,3 hPa (I)
	2,7 hPa (Ex)	3,2 hPa (Ex)
90 l/min	4,0 hPa (I)	5,5 hPa (I)
	4,3 hPa (Ex)	5,4 hPa (Ex)
Moisture loss		
Vt = 500 ml		11 mg/l
Vt = 750 ml		12 mg/l
Vt = 1000 ml		13,9 mg/l

Product features

- Absorb secretions from the tracheostoma quickly to avoid skin irritation
- Act as a pad between the tracheostoma and the tracheostomy tube
- They are slit and can therefore be replaced while wearing the tracheostomy tube
- Available in a variety of materials and absorption strengths:

Nonwoven compresses made of lint-free, skin-friendly material are available in one-layer (A) and two-layer (B) variants (sterile or non sterile)

Compresses coated on one side with aluminium prevent adherence with wound secretion on the skin (C)

Moreover, high-quality tracheal compresses made of special polyurethane foam (PU) are also available which are particularly well tolerated by the skin and can absorb 20 times their own weight in liquid (D, E, F).

These PU compresses are available in a variety of designs and sizes, including in smaller format for children (F)





Indication

- Tracheostomy/ Laryngectomy

		PU	Article number
A	Tracheal compress, one-layer, 9 x 10 cm (slit, with push-through opening), individually packaged, sterile	10 pieces	200230
A	Tracheal compress, one-layer, 9 x 10 cm (slit, with push-through opening), non sterile	10 pieces	200212
B	Tracheal compress, two-layer, 9 x 10 cm (slit, with push-through opening), individually packaged, sterile	10 pieces	200231
B	Tracheal compress, two-layer, 9 x 10 cm (slit, with push-through opening), non sterile	10 pieces	200213
C	Tracheal compress, aluminium-coated, 8 x 9 cm (slit, with push-through opening), individually packaged, sterile	10 pieces	200233
D	Special compress, PU-foam, (slit. with push-through opening), individually packaged, sterile	10 pieces	21858
E	Tracheal compress, PU-foam, 9,6 x 10 cm (slit. with push-through opening), individually packaged, sterile	10 pieces	21840
F	Tracheal compress for children, PU-foam, 6,5 x 6,5 cm, (slit. with push-through opening), individually packaged, sterile	10 pieces	200234

Product features














- For securely holding a tracheostomy tube, a laryngectomy tube or a Stoma button with eyelets in place
- Provide an optimum wearing comfort, as they are made of soft, skin-compatible material (soft fabric, foam, Polyurethan (PU))
- Available in different widths and with different fasteners (hook and slide (narrow and wide), velcro fastener, velcro/ hook)
- Individually adjustable
- Variants for adults and children
- PRIMEDIFIX® with hook enables fast fixation to the eyelets of a tube flange or Stoma Button (A-E)
- PRIMEDIFIX® with velcro fastener is especially suitable for sensitive skin (F-I)



- EASYFIX® with Velcro fastener on one side and a hook on the other enables even motor impaired patients (e.g. after a stroke) to easily fix the neck strap to a tube or stoma button, even with just one hand (J-M)

Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI

		PU	Article number
A	PRIMEDIFIX® neck strap, 1-piece, white, with narrow hook, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200229
B	PRIMEDIFIX® neck strap, 1-piece, skin-colored, with narrow hook, width = 12 mm, length= 26-47 cm, non sterile		10 pieces 200227
C	PRIMEDIFIX® neck strap, 1-piece, transparent, with wide hook, Polyurethan (PU), width = 10 mm, length = 27-50 cm, non sterile		10 pieces 200219
D	PRIMEDIFIX® neck strap, 1-piece, skin-colored, with wide hook, width = 12 mm, length= 26-47 cm, non sterile		10 pieces 200226
E	PRIMEDIFIX® neck strap, 1-piece, white, with wide hook, width = 12 mm, length= 26-47 cm, non sterile		10 pieces 200228
F	PRIMEDIFIX® neck strap for children, 2-part, white, with velcro fastener width = 20 mm, length = 24 cm, non sterile		10 pieces 200203
G	PRIMEDIFIX® neck strap, 2-part, white, with velcro fastener width = 24 mm, length = 47 cm, non sterile		10 pieces 200211
H	PRIMEDIFIX® neck strap, 2-part, white, with velcro fastener width = 16 mm, length = 45,5 cm, sterile		10 pieces 21451
I	PRIMEDIFIX® neck strap, 2-part, white, with velcro fastener, width = 26 mm, length = 45,5 cm, sterile		10 pieces 21450
J	EASYFIX® neck strap, 1-piece, skin-colored, with wide hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200216
K	EASYFIX® neck strap, 1-piece, white, with wide hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200217
L	EASYFIX® neck strap, 1-piece, skin-colored, with narrow hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200214
M	EASYFIX® neck strap, 1-piece, white, with narrow hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200215

COUGH CAP

Product features

- Bypasses the patient’s exhaled air while coughing in order to avoid direct expectoration
- Provides a broad opening to ensure a good airflow
- Can easily be attached to a 15 mm connector of a tracheostomy tube

Indication

- Tracheostomy

	PU	Article number
A Cough cap	1 piece	102081

SHOWER CAP

Product features

- Due to downward opening it protects the tracheostomy tube from incoming liquids while showering
- Ensures good air supply
- Allows for attachment to a 15 mm connector of a tracheostomy tube

Indication

- Tracheostomy
(Not to be used for bathing or swimming!)

	PU	Article number
B Shower cap	1 piece	102080

SHOWER GUARD


Product features

- Covers the whole stoma while showering
- Ensures good tissue compatibility (made of soft medical-grade silicone)
- Adjustable length
- Integrated struts stabilize the cover
- With protective edge which fits seamlessly to the skin
- The openings for breathing are located in the lower area of the device which also avoids the ingress of liquids

Indication

- Tracheostomy/ Laryngectomy
(Not to be used for bathing or swimming!)

	PU	Article number
C Shower guard	1 piece	102085



COUGH CAP SHOWER CAP



SHOWER GUARD

SUCTION CATHETER

Product features

- For secretion suctioning in the trachea, mouth and nose
- Made of a particularly soft and flexible material
- With straight atraumatic tip which is ideal for suctioning in the sensitive trachea
- Two lateral eyes at the catheter tip prevent the catheter from adhering to the inside of the trachea

Indication

- Tracheostomy/ Laryngectomy

SUCTION CATHETER		PU	Article number
A	Suction catheter Ch. 10, Length = 250 mm, sterile	100 pieces	20487
B	Suction catheter Ch. 16, Length = 250 mm, sterile	100 pieces	20488

TRACHEAL SUCTION SET

Product features

- Complete set for tracheal suctioning
- Supplied with an extra sealing cap and an extra label for patient information
- Available with or without vacuum control

Indication

- Tracheostomy/ Laryngectomy

TRACHEAL SUCTION SET		PU	Article number
without picture	TRACHEAL SUCTION SET, without vacuum control, sterile	60 pieces	21214
C	TRACHEAL SUCTION SET with vacuum control, sterile	60 pieces	21215

SUCTION CATHETER



TRACHEAL SUCTION SET

PRIMEDICLEAN®

Product features

- Highly effective product for cleaning and disinfecting tracheostomy tubes
- Comes in a 250ml bottle.

Indication

- Suitable for medical-grade polymer and silver tracheostomy tubes

	PU	Article number
PRIMEDICLEAN® (250 ml)	1 bottle	70810



PRIMEDICLEAN®

CLEANING BRUSHES

Product features

- For gentle and effective cleaning of tracheostomy tubes
- Feature an extremely flexible wire netting and have sturdy yet soft nylon bristles
- Wool tip avoids damaging and scratching of the the tube
- Flexible form allows for Individual adjustment
- With plastic grip for safe handling

Indication

- Tracheostomy/ Laryngectomy

Size	suitable for tubes with an inside diameter of	PU	Article number
6	5-7 mm	5 pieces	200220
8	8-9 mm	5 pieces	200221
10	10-11 mm	5 pieces	200222
12	12-13 mm	5 pieces	200223
14	14-15 mm	5 pieces	200224



CLEANING BRUSHES


LUBRICANT GEL

Product features


- OPTILUBE is a sterile lubricant gel which enables an easy insertion of tracheostomy tubes and stoma buttons




	PU	Article number
Optilube (5 g)	50 pieces	70295




Before Laryngectomy/Tracheostomy




After Laryngectomy



After Tracheostomy
(sieved tracheostomy tube with speaking valve)



After Laryngectomy with voice prosthesis
(sieved tracheostomy tube with HME)



After Laryngectomy with voice prosthesis
(Primedistrip® with speaking valve)



After Laryngectomy
(with electronic speech aid)

Pictograms



With speaking option



With low-pressure cuff



With subglottic suction function



Sterilization with ethylene oxide

MANUFACTURER

Primed[®] HALBERSTADT
MEDIZINTECHNIK GMBH
STRASSE DES 20. JULI 1
D-38820 HALBERSTADT
PHONE.: 0049-(0)3941-668-6
FAX: 0049-(0)3941-24565
PRIMED@PRIMED-HALBERSTADT.DE

YOUR DISTRIBUTOR

ALL RIGHTS TO TECHNICAL MODIFICATIONS RESERVED!

REVISION STATUS: 02/2016-09-29