

Mobile donor lounge DL100P

The perfect lounge ...



... for every blood donation center

Features

- **Portable:** The donor lounge is foldable. Easy to carry on rolls.
- **Shock position:** A collapsed donor can be set into a shock position supported by a gas tube.
- **Donor slip:** Easy to clean and disinfect. The material has a net structure and is comfortable for the donor.
- **Head rest:** The upholstery of the lounge is integrated in the unit.
- **Arm rest:** The armrest is adjustable on both sides.



Adjustable
armrest

Shock position
Gas pressure
supported



Technical data

Material:	Steel painted (round pipe) to prevent scratches and accidents
Weight:	16,5 kg
Dimensions:	fold up: 125 x 55 x 22 cm